

BIOGENESIS KABELIAI IR JUTIKLIAI



SPO2 SENSORS



ECG - EKG CABLES




TEMPERATURE PROBES









SpO2 sensors

- ✓ soft TPU cable
- ✓ double shield signal protection
- ✓ ISO 13485, BPF and CE 0537 certificate
- ✓ internal kevlar® protection
- ✓ DECO® available for clip sensors
- ✓ more than 1000 compatible devices
- ✓ 1 year complete warranty
- ✓ product 100% biocompatible
- ✓ látex free



Sensor types

1 year
warranty

- 1 Adult clip 
- 2 Adult silicone 
- 3 Pediatric clip 
- 4 Pediatric silicone 
- 5 Multisite / Neonatal 
- 6 Wrap silicone 

Accessories available for Neonate / Multisite sensors * sold separately


- ★ Silicone strap  **USD 5**
- ★ Ear clip  **USD 5**
- ★ Veterinary clip  **USD 5**

 DECO® compatible




ECG-EKG cables

- ✓ soft TPU cable
- ✓ double shield signal protection
- ✓ ISO 13485, BPF and CE certificates
- ✓ IEC and AHA color configuration available
- ✓ noise reduction internal carbon film
- ✓ more than 700 compatible devices
- ✓ 1 year complete warranty
- ✓ product 100% biocompatible
- ✓ látex free




Ends



warranty
1 year

Cable types


3 or 5 leads only *



3 Leads / leads

5 Leads / leads

10 Leads / leads complete cable only *



Temperature probes

- ✓ soft TPU cable
- ✓ single shield signal protection
- ✓ ISO 13485, BPF and CE 0537 certificates
- ✓ TITANIUM® available for skin probes
- ✓ more than 500 compatible devices
- ✓ 1 year complete warranty
- ✓ product 100% biocompatible
- ✓ látex free

Probe types

Skin TITANIUM® ★ 

Skin TITANIUM® ★ 


Adult esophageal / rectal ★ 

Pediatric esophageal / rectal ★ 


Timpanic ★ 

Air ★ 


Air ★ 

Air ★ 

Air ★ 

Dual Air ★ 

1 year
warranty



Biogenesis[®]
we care about life

Brand	Product	Page
Adib Jatene	One piece EKG 10 lead cable	
Aerotel	One piece EKG 10 lead cable	12
Air Shields / Hill Rom	Temperature probe	
	Temperature probe	
Ami Systems	Complete SpO2 sensor	13
Anamed	Complete SpO2 sensor	
	Short SpO2 sensor / SpO2 intermediary cable / 3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads)	14
	ECG trunk cable / ECG 3 or 5 lead set	15
Atom	Temperature probe	
BCI	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable / 3 or 5 leads one piece ECG cable	16
	ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	17
Bese	Complete SpO2 sensor	
	Short SpO2 sensor / SpO2 intermediary cable	18
Bionet	Complete SpO2 sensor / Short SpO2 sensor	
	SpO2 intermediary cable / 3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	19
	ECG 3 or 5 lead set	
Bosch	One piece EKG 10 lead cable	20
BTL	Complete SpO2 sensor / Short SpO2 sensor	
	SpO2 intermediary cable / One piece EKG 10 lead cable	21
Burdick	One piece EKG 10 lead cable	
Cardioline / Kenz	Cable de ECG de 5 leads / One piece EKG 10 lead cable	22
Cardiotek	One piece EKG 10 lead cable	
Cardiorapid	One piece EKG 10 lead cable	23
Choice	Complete SpO2 sensor	
	Short SpO2 sensor / SpO2 intermediary cable / Temperature probe	24

Brand	Product	Page
Cmos Drake	Complete SpO2 sensor	25
	Short SpO2 sensor / SpO2 intermediary cable	26
Combiomed	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable	27
Corometrics	3 or 5 leads one piece ECG cable	28
	ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	
Criticare CSI	Complete SpO2 sensor	29
	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable	
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	
Critikon	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable	31
Danatech	Complete SpO2 sensor	32
Datascope	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable	
	SpO2 intermediary cable / 3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads)	
Digicare	ECG trunk cable / ECG 3 or 5 lead set	34
	Complete SpO2 sensor	35
Dixtal	Short SpO2 sensor / SpO2 intermediary cable	
	Complete SpO2 sensor / Cable de ECG de 3 o 5 leads	36
	ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set / One piece EKG 10 lead cable	
Dräger / Siemens	Temperature probe	37
	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable	38
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	39
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	40
Ecafix	Temperature probe	41
	Complete SpO2 sensor / 3 or 5 leads one piece ECG cable	42

Brand	Product	Page
Ecafix	ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	43
Edan	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable / One piece EKG 10 lead cable	44
Electronica Trentina / Cardiette	One piece EKG 10 lead cable	45
Emai Transmai	Complete SpO2 sensor / Short SpO2 sensor	46
	SpO2 intermediary cable / 3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	
Esaote	ECG 3 or 5 lead set	47
	3 or 5 leads one piece ECG cable / One piece EKG 10 lead cable	
Fanem	Temperature probe	48
Feas	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	49
	ECG 3 or 5 lead set	
Fischer & Paykel	Temperature probe / Dual temperature probe	50
Fukuda Denshi	Complete SpO2 sensor	51
	Short SpO2 sensor / SpO2 intermediary cable / One piece EKG 10 lead cable	
GE - Critikon	Cable de ECG de 3 o 5 leads	52
	Cable de ECG desmontable de 3 o 5 leads / ECG trunk cable / ECG 3 or 5 lead set	
GE - Datex Ohmeda	Complete SpO2 sensor	53
	Complete SpO2 sensor / Short SpO2 sensor	54
	Short SpO2 sensor / SpO2 intermediary cable	55
	SpO2 intermediary cable / 3 or 5 leads one piece ECG cable	56
	Cable de ECG desmontable de 3 o 5 leads / ECG trunk cable / ECG 3 or 5 lead set	57
GE - Marquette	Temperature probe	58
	Complete SpO2 Sensor	59
	Complete SpO2 sensor / Short SpO2 sensor	60
	SpO2 intermediary cable / 3 or 5 leads one piece ECG cable	61
	Cable de ECG desmontable de 3 o 5 leads / ECG trunk cable / ECG 3 or 5 lead set	62
GE - Marquette	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	62

Brand	Product	Page
GE - Marquette	One piece EKG 10 lead cable / Temperature probe	63
General Meditech	Complete SpO2 sensor	64
	Complete SpO2 sensor	
Innomed	One piece EKG 10 lead cable	64
Instramed	Complete SpO2 sensor	65
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	
Invivo Research	Complete SpO2 sensor / Short SpO2 sensor	66
	SpO2 intermediary cable	67
Ivy Biomedical	Complete SpO2 sensor	
Kenz	One piece EKG 10 lead cable	68
Kontron	Complete SpO2 sensor	
	Short SpO2 sensor / SpO2 intermediary cable / One piece EKG 10 lead cable	
Ktmed Epox	Complete SpO2 sensor	69
	Short SpO2 sensor / SpO2 intermediary cable	
MDE	Complete SpO2 sensor	70
Medexel	Complete SpO2 sensor	
Medix	Temperature probe	71
Medinova	Complete SpO2 sensor	
Medlab	Complete SpO2 sensor	71
	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable	72
	SpO2 intermediary cable / 3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	73
	ECG 3 or 5 lead set	74
Mennen	Complete SpO2 sensor / Short SpO2 sensor	75
	SpO2 intermediary cable / 3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads)	
	ECG trunk cable / ECG 3 or 5 lead set / 3 or 5 leads one piece ECG cable	76

Brand	Product	Page
Mennen	ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set / Temperature probe	77
Mindray	Complete SpO2 sensor	78
	Short SpO2 sensor / SpO2 intermediary cable	79
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	80
	ECG 3 or 5 lead set / Temperature probe	81
Moriya	Complete SpO2 sensor / Short SpO2 sensor / 3 or 5 leads one piece ECG cable	82
	ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	83
Mortara	One piece EKG 10 lead cable	84
	One piece EKG 10 lead cable	84
Nellcor	Complete SpO2 sensor	85
	Complete SpO2 sensor	85
	Short SpO2 sensor / SpO2 intermediary cable	86
	SpO2 intermediary cable	87
Nihon Kohden	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	88
	Complete SpO2 sensor / Short SpO2 sensor	89
	SpO2 intermediary cable	90
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	91
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	92
Nonin	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	93
	ECG 3 or 5 lead set / One piece EKG 10 lead cable	94
	Temperature probe	95
Novamatrix	Complete SpO2 sensor / Short SpO2 sensor	96
	SpO2 intermediary cable	96
Novamatrix	Complete SpO2 sensor	96
	Short SpO2 sensor / SpO2 intermediary cable	97

Brand	Product	Page
Ohio	Temperature probe	97
Pace Tech	Complete SpO2 sensor	98
	Short SpO2 sensor / SpO2 intermediary cable	99
Palco	Complete SpO2 sensor	100
Philips - HP	Complete SpO2 sensor	100
	Short SpO2 sensor / SpO2 intermediary cable	101
	SpO2 intermediary cable / 3 or 5 leads one piece ECG cable / ECG trunk cable de 3 y 5 leads	102
	ECG trunk cable / ECG 3 or 5 lead set / 3 or 5 leads one piece ECG cable	103
	ECG trunk cable (3 or 5 leads) / ECG trunk cable	104
	ECG 3 or 5 lead set / One piece EKG 10 lead cable / Temperature probe	105
Physio Control / Medtronic	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	106
	ECG 3 or 5 lead set / 3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads)	107
	ECG trunk cable / ECG 3 or 5 lead set	108
	One piece EKG 10 lead cable	108
Progetti	Complete SpO2 sensor	109
	Short SpO2 sensor / SpO2 intermediary cable	109
Rosch	Temperature probe	109
Schiller	One piece EKG 10 lead cable	110
	One piece EKG 10 lead cable	110
Spacelabs	Complete SpO2 sensor	111
	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable	111
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable	112
	ECG 3 or 5 lead set / 3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads)	113
	ECG trunk cable / ECG 3 or 5 lead set	114

Brand	Product	Page
Takaoka	Complete SpO2 sensor	114
	Short SpO2 sensor / SpO2 intermediary cable	115
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	116
TEB	Complete SpO2 sensor / Short SpO2 sensor / SpO2 intermediary cable	117
	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads) / ECG trunk cable / ECG 3 or 5 lead set	118
Welch Allyn	Temperature probe / 3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads)	119
	ECG trunk cable / ECG 3 or 5 lead set	120
YSI 400	Temperature probe	120
YSI 700	Temperature probe	121
Zoll	3 or 5 leads one piece ECG cable / ECG trunk cable (3 or 5 leads)	121
	ECG trunk cable / ECG 3 or 5 lead set	122

Important


All BIOGENESIS products here listed are compatible with the mentioned brands, and those brands are trademarks of their respective companies.


This products are manufactured in the Montevideo Biogénesis facilities under the supervision of expert engineers and tested under the most stringent quality standards in compliance with the ISO 9001, ISO 13485, CE, CE 0537, and good manufacturing practices.

Biogénesis S.R.L. team.


Adib Jatene

One piece EKG
10 lead cable.

Options		Code
Plug / Banana		BIOECG.AJ.P10F




< Cable length: 3,0m >




Compatible with: All Adib Jatene devices.


Aerotel

One piece EKG
10 lead cable.

Options		Code
Plug / Banana		BIOECG.AE.P10F





< Cable length: 3,0m >



Compatible with: All Aerotel devices.


Air Shields / Hill Rom

Temperature probe.

Options		Code
Skin TITANIUM® 9,5mm		BIOTEMP.200PT
Skin TITANIUM® 5,5mm		BIOTEMP.200PTP




< Cable length: 1,5m >



Compatible with: C-100, C-100QT, C-200, C-450, C-450QT, 550QT; Warmer Systems: 7810, 7820, 7830, 7840.

Air Shields / Hill Rom

Temperature probe.

Options		Code
Skin TITANIUM® 9,5mm		BIOTEMP.700PT




< Cable length: 1,5m >



Compatible with: Isolette C2000, Versalet 7700.


Air Shields / Hill Rom

Temperature probe.

Options		Code	
Air		BIOTEMP.200A4	




< Cable length: 1,5m >



Compatible with: C450.


Air Shields / Hill Rom

Temperature probe.

Options		Code	
Air		BIOTEMP.250A4	



< Cable length: 1,5m >



Compatible with: C-100, C-100QT, C-200, C-450.

AMI Systems







SpO2 sensor.

Compatible with:
AMI 7300 Series.



< Cable length: 2,7m >



Options		Code	
Adult Clip		BIOOX.AM.AC3	
Adult Silicone		BIOOX.AM.AS3	
Pediatric Clip		BIOOX.AM.PC3	
Pediatric Silicone		BIOOX.AM.PS3	
Neonate / Multisite		BIOOX.AM.Y3	
Wrap silicone		BIOOX.AM.WS3	


Anamed







SpO2 sensor.


Compatible with:
78100-B, 78600, Vital Line, Vital Wave.




< Cable length: 2,7m >



Options		Code	
Adult Clip		BIOOX.PR.AC3	
Adult Silicone		BIOOX.PR.AS3	
Pediatric Clip		BIOOX.PR.PC3	
Pediatric Silicone		BIOOX.PR.PS3	
Neonate / Multisite		BIOOX.PR.Y3	
Wrap silicone		BIOOX.PR.WS3	



< Cable length: 0,9m >




Anamed

SpO2 short sensor.
this product requires an extension cable

Compatible with:
78100-B, 78600, Vital Line, Vital Wave.

Options		Code	
Adult Clip		BIOOX.BC.AC1	
Adult Silicone		BIOOX.BC.AS1	
Pediatric Clip		BIOOX.BC.PC1	
Pediatric Silicone		BIOOX.BC.PS1	
Neonate / Multisite		BIOOX.BC.Y1	
Wrap silicone		BIOOX.BC.WS1	

Anamed



< Cable length: 2,0m >



SpO2 intermediary cable.

Compatible with:
78100-B, 78600, Vital Line, Vital Wave.


Description	Code	
SpO2 intermediary cable	BIOOX.PR.IC2	

Anamed


One piece ECG cable (3 or 5 leads)

Compatible with:
AM-78100A, AM 78100B,
AM 78600, Vital Lite,
Vital Wave, AM-7810A.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



< Cable length: 3,5m >




Anamed


ECG trunk cable with leads (3 or 5 leads)

Compatible with:
AM-78100A, AM 78100B,
AM 78600, Vital Lite,
Vital Wave, AM-7810A.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




Anamed

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	







< Cable length: 3,0m >



Compatible with: AM-78100A, AM 78100B, AM 78600, Vital Lite, Vital Wave, AM-7810A.


Anamed

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	



Compatible with:
AM-78100A, AM 78100B,
AM 78600, Vital Lite,
Vital Wave, AM-7810A.

Cable length: 1,0m >




Atom

Temperature probe.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.500PT	
Skin TITANIUM® 5,5mm		BIOTEMP.500PTP	





< Cable length: 1,5m >



Compatible with: All Atom devices with Plug 1/4.

Atom

Temperature probe.


Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.550PT	
Skin TITANIUM® 5,5mm		BIOTEMP.550PTP	




< Cable length: 1,5m >



Compatible with: V-2100, 2100G, 2200, 3200, 5005, 505ST, 80, 800, 808, 850 incubadoras o infant warmers.



< Cable length: 2,7m >



BCI

SpO2 sensor.

Compatible with:

3300, 3301, 3302, 3303, 3304, 3400, 3401, 3403, 9200; Spacelabs Burdick: PRO2; Drager: Oxipac 2000; MedChoice: MD300A; Invacare 3300nv.

Options		Code	
Adult Clip		BIOOX.BC.AC3	
Adult Silicone		BIOOX.BC.AS3	
Pediatric Clip		BIOOX.BC.PC3	
Pediatric Silicone		BIOOX.BC.PS3	
Neonate / Multisite		BIOOX.BC.Y3	
Wrap silicone		BIOOX.BC.WS3	



< Cable length: 0,9m >



BCI

SpO2 short sensor.


this product requires an extension cable

Compatible with:

3180, 3401, 3304, 3303, 3302, 3301, 3300, 3400, 3403, Smiths Medical BCI Capnocheck 9004, Spectro 2 | 10, 20, 30, SurgiVet; Spacelabs Burdick: PRO2; Drager: Drager Oxipac 2000; Armstrong Medical AD-1000; Beijing Choice Electronic: Pulse Oximeter model MD30011.

Options		Code	
Adult Clip		BIOOX.BC.AC1	
Adult Silicone		BIOOX.BC.AS1	
Pediatric Clip		BIOOX.BC.PC1	
Pediatric Silicone		BIOOX.BC.PS1	
Neonate / Multisite		BIOOX.BC.Y1	
Wrap silicone		BIOOX.BC.WS1	

BCI



sensor connection

< Cable length: 2,0m >



device connection

SpO2 intermediary cable.

Compatible with:

All BCI devices.

Description	Code	
SpO2 intermediary cable	BIOOX.BC.IC2	


BCI

One piece ECG cable (3 or 5 leads)


Compatible with:

3101, 3404, 6100, 6200, 9100, 9200, Advisor, Autocar, Life Span 100 (2005).

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



< Cable length: 3,5m >







BCI

ECG trunk cable with leads
(3 or 5 leads)


Compatible with:

3101, 3404, 6100, 6200, 9100,
9200, Advisor, Autocar, Life
Span 100 (2005).

Options		Code
3 leads snap		BIOECG.BC.S3D
3 leads clamp		BIOECG.BC.C3D
5 leads snap		BIOECG.BC.S5D
5 leads clamp		BIOECG.BC.C5D






< Cable length: 4,0m >




BCI

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code
3 leads		BIOECG.BC.3T
5 leads		BIOECG.BC.5T



< Cable length: 3,0m >



Compatible with:





3101, 3404, 6100, 6200, 9100, 9200, Advisor, Autocar, Life Span 100 (2005).

BCI


ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:


3101, 3404, 6100, 6200, 9100,
9200, Advisor, Autocar, Life
Span 100 (2005).

Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5

Cable length: 1,0m >



< Cable length: 2,7m >










Bese

SpO2 sensor.


Compatible with:

Biomonitor 4, Biomonitor 7.

Options		Code
Adult Clip		BIOOX.BE.AC3
Adult Silicone		BIOOX.BE.AS3
Pediatric Clip		BIOOX.BE.PC3
Pediatric Silicone		BIOOX.BE.PS3
Neonate / Multisite		BIOOX.BE.Y3
Wrap silicone		BIOOX.BE.WS3



< Cable length: 0,9m >









Bese

SpO2 short sensor.


this product requires an extension cable

Compatible with:

Biomonitor 4, Biomonitor 7.


Options		Code	
Adult Clip		BIOOX.BC.AC1	
Adult Silicone		BIOOX.BC.AS1	
Pediatric Clip		BIOOX.BC.PC1	
Pediatric Silicone		BIOOX.BC.PS1	
Neonate / Multisite		BIOOX.BC.Y1	
Wrap silicone		BIOOX.BC.WS1	

Bese



sensor connection

< Cable length: 2,0m >



device connection

SpO2 intermediary cable.

Compatible with:

Biomonitor 4, Biomonitor 7.

Description	Code	
SpO2 intermediary cable	BIOOX.BE.IC2	



< Cable length: 2,7m >









Bionet

SpO2 sensor.

Compatible with:

BM1, BM3, BM3 Plus.

Options		Code	
Adult Clip		BIOOX.BI.AC3	
Adult Silicone		BIOOX.BI.AS3	
Pediatric Clip		BIOOX.BI.PC3	
Pediatric Silicone		BIOOX.BI.PS3	
Neonate / Multisite		BIOOX.BI.Y3	
Wrap silicone		BIOOX.BI.WS3	



< Cable length: 0,9m >









Bionet


SpO2 short sensor.

this product requires an extension cable

Compatible with:

BM1, BM3, BM3 Plus.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Bionet

SpO2 intermediary cable.

Compatible with:

BM1, BM3, BM3 Plus.





Description	Code	
SpO2 intermediary cable	BIOOX.BE.IC2	


Bionet

One piece ECG cable
(3 or 5 leads)


Compatible with:

BM3.

Options		Code	
3 leads snap		BIOECG.BI.S3F	
3 leads clamp		BIOECG.BI.C3F	
5 leads snap		BIOECG.BI.S5F	
5 leads clamp		BIOECG.BI.C5F	



< Cable length: 3,5m >








Bionet

ECG trunk cable with leads
(3 or 5 leads)

Compatible with:

BM3.

Options		Code	
3 leads snap		BIOECG.BI.S3D	
3 leads clamp		BIOECG.BI.C3D	
5 leads snap		BIOECG.BI.S5D	
5 leads clamp		BIOECG.BI.C5D	





< Cable length: 4,0m >




Bionet

ECG trunk cable
(3 or 5 leads).

this product requires a lead set

Options		Code	
3 leads		BIOECG.BI.3T	
5 leads		BIOECG.BI.5T	



< Cable length: 3,0m >







Compatible with:

BM3.

Bionet

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
BM3.


Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >




Bosch


One piece EKG
10 lead cable.




< Cable length: 3,0m >



Compatible with: 603 D, 603 ER, 603 S, 606.









< Cable length: 2,7m >




BTL


SpO2 sensor.

Compatible with:
All BTL devices.

Options		Code	
Adult Clip		BIOOX.BL.AC3	
Adult Silicone		BIOOX.BL.AS3	
Pediatric Clip		BIOOX.BL.PC3	
Pediatric Silicone		BIOOX.BL.PS3	
Neonate / Multisite		BIOOX.BL.Y3	
Wrap silicone		BIOOX.BL.WS3	









< Cable length: 0,9m >




BTL

SpO2 short sensor.
this product requires an extension cable


Compatible with:
All BTL devices.

Options		Code	
Adult Clip		BIOOX.BLD.AC1	
Adult Silicone		BIOOX.BLD.AS1	
Pediatric Clip		BIOOX.BLD.PC1	
Pediatric Silicone		BIOOX.BLD.PS1	
Neonate / Multisite		BIOOX.BLD.Y1	
Wrap silicone		BIOOX.BLD.WS1	



sensor connection

< Cable length: 2,0m >



device connection

BTL


SpO2 intermediary cable.

Compatible with:
All BTL devices.


Description	Code	
SpO2 intermediary cable	BIOOX.BL.IC2	

BTL


One piece EKG
10 lead cable.



< Cable length: 3,0m >




Compatible with: All BTL devices.


Options	Code	
Plug / Banana 	BIOECG.SC.P10F	

Burdick


One piece EKG
10 lead cable.



< Cable length: 3,0m >




Compatible with: EK10, Elite, E350, E350I, E550, E560, Eclipse LE, LE I, Elite II, Quinton Eclipse Plus.

Options	Code	
Plug / Banana 	BIOECG.BU.P10F	

Burdick


One piece EKG
10 lead cable.



< Cable length: 3,0m >






Compatible with: Eclipse 4, 300, 400, 800, 850, Plus and Premier, Elite II, Eclipse LE II and Atria Systems, Atria 3000, 3100, 6100.

Options	Code	
Plug / Banana 	BIOECG.BU2.P10F	

Cardioline

One piece ECG cable
5 leads only.

Options		Code	
5 leads snap		BIOECG.CDK.S5F	
5 leads clamp		BIOECG.SDK.S5F	




< Cable length: 0,6m >



Compatible with: Click Holter.


Cardioline

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.CAK.P10F	


Compatible with:
Click Holter.


Cable length: 0,6m >




Cardioline

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.NK.P10F	




< Cable length: 3,0m >



Compatible with: E1, ER1, E3, ITA300, C3.


Cardioline / Kenz

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.CA.P10F	




< Cable length: 3,0m >




Compatible with: ITA 1 -TA1 Plus, ITA 3 - ITA 3 Plus, ITA 30, ITA 60 - ITA 60 Plus, Elan, Elan 100, Elan 1100, 103, 106, 107, 302, 303, 1201 and Cardioline Remco Delta 1, Delta 3, Delta 30, Delta 30D, Delta 60, Delta 1 Plus, Delta 3 Plus, Delta 60 Plus, ECG 106.


Cardiotek

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.SC.P10F	




< Cable length: 3,0m >




Compatible with: EP-Tracer series.


Cardiorapid

One piece EKG
10 lead cable.


Options		Code	
Plug / Banana		BIOECG.SC.P10F	




< Cable length: 3,0m >



Compatible with: K11, K111, K112, K131, K31, K33, K36.









< Cable length: 2,7m >




Choice

SpO2 sensor.

Compatible with:
All Choice
devices with Redel 5 pin connector.

Options		Code	
Adult Clip		BIOOX.CH.AC3	
Adult Silicone		BIOOX.CH.AS3	
Pediatric Clip		BIOOX.CH.PC3	
Pediatric Silicone		BIOOX.CH.PS3	
Neonate / Multisite		BIOOX.CH.Y3	
Wrap silicone		BIOOX.CH.WS3	









< Cable length: 2,7m >




Choice


SpO2 sensor.

Compatible with:
MD300A.

Options		Code	
Adult Clip		BIOOX.CH2.AC3	
Adult Silicone		BIOOX.CH2.AS3	
Pediatric Clip		BIOOX.CH2.PC3	
Pediatric Silicone		BIOOX.CH2.PS3	
Neonate / Multisite		BIOOX.CH2.Y3	
Wrap silicone		BIOOX.CH2.WS3	



< Cable length: 0,9m >




Choice

SpO2 short sensor.
this product requires an extension cable


Compatible with:
All Choice devices.

Options		Code	
Adult Clip		BIOOX.BC.AC1	
Adult Silicone		BIOOX.BC.AS1	
Pediatric Clip		BIOOX.BC.PC1	
Pediatric Silicone		BIOOX.BC.PS1	
Neonate / Multisite		BIOOX.BC.Y1	
Wrap silicone		BIOOX.BC.WS1	

Choice



< Cable length: 2,0m >




SpO2 intermediary cable.

Compatible with:
All Choice
devices with Redel 5 pin connector.

Description	Code	
SpO2 intermediary cable	BIOOX.CH.IC2	

Choice



< Cable length: 2,0m >



SpO2 intermediary cable.

Compatible with:
MD300A.


Description	Code	
SpO2 intermediary cable	BIOOX.CH2.IC2	

Choice

Temperature probe.



< Cable length: 3,0m >



Compatible with:
All Choice devices.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.400PT	
Skin TITANIUM® 5,5mm		BIOTEMP.400PTP	
Timpanic		BIOTEMP.400T	
Adult esophageal / rectal		BIOTEMP.400RA	
Pediatric esophageal / rectal		BIOTEMP.400RP	









< Cable length: 2,7m >



Cmos Drake

SpO2 sensor.

Compatible with:
All Cmos Drake
devices with Redel 5 pin connector.

Options		Code	
Adult Clip		BIOOX.CH.AC3	
Adult Silicone		BIOOX.CH.AS3	
Pediatric Clip		BIOOX.CH.PC3	
Pediatric Silicone		BIOOX.CH.PS3	
Neonate / Multisite		BIOOX.CH.Y3	
Wrap silicone		BIOOX.CH.WS3	









< Cable length: 2,7m >



Cmos Drake


SpO2 sensor.

Compatible with:
Oxiline Plus.

Options		Code	
Adult Clip		BIOOX.CMO.AC3	
Adult Silicone		BIOOX.CMO.AS3	
Pediatric Clip		BIOOX.CMO.PC3	
Pediatric Silicone		BIOOX.CMO.PS3	
Neonate / Multisite		BIOOX.CMO.Y3	
Wrap silicone		BIOOX.CMO.WS3	









< Cable length: 2,7m >



Cmos Drake


SpO2 sensor.

Compatible with:
All Cmos Drake
devices with DB9 connector.

Options		Code	
Adult Clip		BIOOX.BC.AC3	
Adult Silicone		BIOOX.BC.AS3	
Pediatric Clip		BIOOX.BC.PC3	
Pediatric Silicone		BIOOX.BC.PS3	
Neonate / Multisite		BIOOX.BC.Y3	
Wrap silicone		BIOOX.BC.WS3	









< Cable length: 0,9m >




Cmos Drake

SpO2 short sensor.
this product requires an extension cable

Compatible with:
All Cmos Drake
devices with DB9 connector.

Options		Code	
Adult Clip		BIOOX.BC.AC1	
Adult Silicone		BIOOX.BC.AS1	
Pediatric Clip		BIOOX.BC.PC1	
Pediatric Silicone		BIOOX.BC.PS1	
Neonate / Multisite		BIOOX.BC.Y1	
Wrap silicone		BIOOX.BC.WS1	

Cmos Drake



sensor connection

< Cable length: 2,0m >




device connection

SpO2 intermediary cable.

Compatible with:
All Cmos Drake
devices with Redel 5 pin connector.

Description	Code	
SpO2 intermediary cable	BIOOX.CH.IC2	

Cmos Drake



sensor connection

< Cable length: 2,0m >




device connection

SpO2 intermediary cable.

Compatible with:
Oxilife Plus.

Description	Code	
SpO2 intermediary cable	BIOOX.CMO.IC2	

Cmos Drake



sensor connection

< Cable length: 2,0m >



device connection

SpO2 intermediary cable.

Compatible with:
All Cmos Drake
devices with DB9 connector.

Description	Code	
SpO2 intermediary cable	BIOOX.BC.IC2	



< Cable length: 2,7m >




Combiomed


SpO2 sensor.

Compatible with:
All Combiomed devices.

Options		Code	
Adult Clip		BIOOX.NN.AC3	
Adult Silicone		BIOOX.NN.AS3	
Pediatric Clip		BIOOX.NN.PC3	
Pediatric Silicone		BIOOX.NN.PS3	
Neonate / Multisite		BIOOX.NN.Y3	
Wrap silicone		BIOOX.NN.WS3	



< Cable length: 0,9m >




Combiomed

SpO2 short sensor.
this product requires an extension cable

Compatible with:
All Combiomed devices.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Combiomed

SpO2 intermediary cable.

Compatible with: All Combiomed devices.


Description	Code	
SpO2 intermediary cable	BIOOX.NN.IC2	

Corometrics


One piece ECG cable
(3 or 5 leads)

Compatible with:
556.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	







< Cable length: 3,5m >




Corometrics


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
556.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




Corometrics

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	



< Cable length: 3,0m >







Compatible with: 556.

Corometrics

ECG lead set
(3 or 5 leads)
this product requires a trunk cable


Compatible with:
556.

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >









< Cable length: 2,7m >




Criticare CSI


SpO2 sensor.

Compatible with:
Scholar, Poet (602-3A), 1100, 2200, 500, 500A
501, 502 series, 503, 504, 504/P, 504O/P,
504US/P, 507, 507E/P/R, 507/S, 507SD, 507O/P,
508, 509, 581.

Options		Code	
Adult Clip		BIOOX.CR.AC3	
Adult Silicone		BIOOX.CR.AS3	
Pediatric Clip		BIOOX.CR.PC3	
Pediatric Silicone		BIOOX.CR.PS3	
Neonate / Multisite		BIOOX.CR.Y3	
Wrap silicone		BIOOX.CR.WS3	



< Cable length: 2,7m >



Criticare CSI

SpO2 sensor.


Compatible with:

Alaris 4415 (IVAC): VitalCheck 4400, 4410; Criticare: 503DX-W, 503S, 503DX, 503 Spot, 504DX/P, 504DX, 506 series (except 506LN3), 507ELC2, 8000 series, 8100 series, Poet Plus, Poet R Plus, Poet (602- 3A), Scholar, Scholar III, nGenuity.

Options		Code	
Adult Clip		BIOOX.CRD.AC3	
Adult Silicone		BIOOX.CRD.AS3	
Pediatric Clip		BIOOX.CRD.PC3	
Pediatric Silicone		BIOOX.CRD.PS3	
Neonate / Multisite		BIOOX.CRD.Y3	
Wrap silicone		BIOOX.CRD.WS3	



< Cable length: 0,9m >



Criticare CSI

SpO2 short sensor.


this product requires an extension cable

Compatible with:

Alaris (IVAC): VitalCheck 4400 511SD; Criticare: 503DX-W, 503S, 503DX, 503 Spot, Mini Spot, 504DX/P, 504DX, 506 series (except 506LN3), 507ELC2, 507EL, 8100 series, Poet Plus, Poet (602- 3A), Scholar, Scholar III.

Options		Code	
Adult Clip		BIOOX.CRD.AC1	
Adult Silicone		BIOOX.CRD.AS1	
Pediatric Clip		BIOOX.CRD.PC1	
Pediatric Silicone		BIOOX.CRD.PS1	
Neonate / Multisite		BIOOX.CRD.Y1	
Wrap silicone		BIOOX.CRD.WS1	

Criticare CSI



sensor connection

< Cable length: 2,0m >




device connection

SpO2 intermediary cable.

Compatible with:


Scholar, Poet (602-3A), 1100, 581, 508, 507, 507SD, 504, 503, 501; Mindray: 500A.

Description	Code	
SpO2 intermediary cable	BIOOX.CR.IC2	



sensor connection

< Cable length: 2,0m >



device connection

Criticare CSI

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.CRD.IC2	





Compatible with: 506 Spot CSI, 507SD, 503 Spot.


Criticare CSI

One piece ECG cable (3 or 5 leads)


Compatible with:

504US, 506, 507, 508, 602-4USP, 1100, 1508, 2200, Poet 2, Scholar.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



< Cable length: 3,5m >








Criticare CSI

ECG trunk cable with leads (3 or 5 leads)


Compatible with:

504US, 506, 507, 508, 602-4USP, 1100, 1508, 2200, Poet 2, Scholar.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	





< Cable length: 4,0m >




Criticare CSI


ECG trunk cable (3 or 5 leads).

this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	



< Cable length: 3,0m >



Compatible with: 504US, 506, 507, 508, 602-4USP, 1100, 1508, 2200, Poet 2, Scholar.

Criticare CSI





ECG lead set (3 or 5 leads)

this product requires a trunk cable

Compatible with:

504US, 506, 507, 508, 602-4USP, 1100, 1508, 2200,

Poet 2, Scholar.

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >





< Cable length: 2,7m >









Critikon

SpO2 sensor.

Compatible with:

Dinamap 8700, 8710, 8720 Pulsse, 8721, 9710, 9700, 8280 and Oxyshuttle.

Options		Code	
Adult Clip		BIOOX.CK.AC3	
Adult Silicone		BIOOX.CK.AS3	
Pediatric Clip		BIOOX.CK.PC3	
Pediatric Silicone		BIOOX.CK.PS3	
Neonate / Multisite		BIOOX.CK.Y3	
Wrap silicone		BIOOX.CK.WS3	



< Cable length: 0,9m >









Critikon

SpO2 short sensor.


this product requires an extension cable

Compatible with:

Dinamap 9700, 9710, 9720 Plus, Compact, MPS, Select, Dinamap Pro 100-400.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	

Critikon



sensor connection

< Cable length: 2,0m >



device connection

SpO2 intermediary cable.

Compatible with:

Dinamap 8700, 8710, 8720 Pulsse, 8721, 9710, 9700, 8280 and Oxyshuttle.

Description	Code	
SpO2 intermediary cable	BIOOX.CK.IC2	



< Cable length: 2,7m >









Danatech

SpO2 sensor.

Compatible with:

All Danatech devices.

Options		Code	
Adult Clip		BIOOX.MI.AC3	
Adult Silicone		BIOOX.MI.AS3	
Pediatric Clip		BIOOX.MI.PC3	
Pediatric Silicone		BIOOX.MI.PS3	
Neonate / Multisite		BIOOX.MI.Y3	
Wrap silicone		BIOOX.MI.WS3	



< Cable length: 2,7m >



Datascope

SpO2 sensor.

Compatible with:

All Datascope devices with Redel 6 pin connector.

Options		Code	
Adult Clip		BIOOX.MI.AC3	
Adult Silicone		BIOOX.MI.AS3	
Pediatric Clip		BIOOX.MI.PC3	
Pediatric Silicone		BIOOX.MI.PS3	
Neonate / Multisite		BIOOX.MI.Y3	
Wrap silicone		BIOOX.MI.WS3	



< Cable length: 2,7m >




Datascope

SpO2 sensor.


Compatible with:

Accusat, Accutorr 3 SAT, Accutorr 4 SAT, Accutorr Plus, Passport, Passport 2, Passport 2LT, Passport 5L, Passport IL, Passport EL, Passport XL, Passport XG, Datascope Multinex Plus Monitor, Datascope 3000, Passport D.

Options		Code	
Adult Clip		BIOOX.DS.AC3	
Adult Silicone		BIOOX.DS.AS3	
Pediatric Clip		BIOOX.DS.PC3	
Pediatric Silicone		BIOOX.DS.PS3	
Neonate / Multisite		BIOOX.DS.Y3	
Wrap silicone		BIOOX.DS.WS3	



< Cable length: 0,9m >



Dastascope

SpO2 short sensor.


this product requires an extension cable

Compatible with:

Accusat, Accutorr 3 SAT, Accutorr 4 SAT, Accutorr Plus, Passport, Passport 2, Passport 2LT, Passport IL, Passport XL and Passport XG.


Options		Code	
Adult Clip		BIOOX.DSD.AC1	
Adult Silicone		BIOOX.DSD.AS1	
Pediatric Clip		BIOOX.DSD.PC1	
Pediatric Silicone		BIOOX.DSD.PS1	
Neonate / Multisite		BIOOX.DSD.Y1	
Wrap silicone		BIOOX.DSD.WS1	

Datascope



sensor connection

< Cable length: 2,0m >




device connection

SpO2 intermediary cable.

Compatible with:

All Datascope devices with Redel 6 pin connector.

Description	Code	
SpO2 intermediary cable	BIOOX.MI.IC2	



Datascope

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.DS.IC2	


Compatible with:

Accusat, Accutorr 3 SAT, Accutorr 4 SAT, Accutorr Plus, Passport, Passport 2, Passport 2LT, Passport IL, Passport XL and Passport XG.

Datascope

One piece ECG cable
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



Compatible with:

Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Ballon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96, System 97.

Datascope

ECG trunk cable with leads
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






Compatible with:

Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Ballon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96, System 97.


Datascope

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code
3 leads		BIOECG.BC.3T
5 leads		BIOECG.BC.5T







< Cable length: 3,0m >



Compatible with: Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Ballon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96, System 97.

Datascope

ECG lead set
(3 or 5 leads)
this product requires a trunk cable


Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5

Cable length: 1,0m >



Compatible with:

Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Ballon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96, System 97.









< Cable length: 2,7m >



Digicare

SpO2 sensor.

Compatible with:
PO 900, PO 920, PO 930, 5000, DIGIPRESS.

Options		Code
Adult Clip		BIOOX.DG.AC3
Adult Silicone		BIOOX.DG.AS3
Pediatric Clip		BIOOX.DG.PC3
Pediatric Silicone		BIOOX.DG.PS3
Neonate / Multisite		BIOOX.DG.Y3
Wrap silicone		BIOOX.DG.WS3



< Cable length: 0,9m >




Digicare

SpO2 short sensor.
this product requires an extension cable


Compatible with:
PO 900, PO 920, PO 930, 5000, DIGIPRESS.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Digicare

SpO2 intermediary cable.

Compatible with:
PO 900, PO 920, PO 930, 5000, DIGIPRESS.

Description	Code	
SpO2 intermediary cable	BIOOX.DG.IC2	



< Cable length: 2,7m >



Dixtal

SpO2 sensor.

Compatible with:
DX 405, DX 2010, DX 2020, DX 2021,
DX 2022, DX 2023, DX 2405, DX 2515,
DX 7100, DX 7300, DX 8100.


Options		Code	
Adult Clip		BIOOX.DI.AC3	
Adult Silicone		BIOOX.DI.AS3	
Pediatric Clip		BIOOX.DI.PC3	
Pediatric Silicone		BIOOX.DI.PS3	
Neonate / Multisite		BIOOX.DI.Y3	
Wrap silicone		BIOOX.DI.WS3	

Dixtal


One piece ECG cable
(3 or 5 leads)

Compatible with:
All Dixtal devices.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	







< Cable length: 3,5m >




Dixtal


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
All Dixtal devices.

Options		Code
3 leads snap		BIOECG.BC.S3D
3 leads clamp		BIOECG.BC.C3D
5 leads snap		BIOECG.BC.S5D
5 leads clamp		BIOECG.BC.C5D






< Cable length: 4,0m >




Dixtal

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code
3 leads		BIOECG.BC.3T
5 leads		BIOECG.BC.5T



< Cable length: 3,0m >







Compatible with: All Dixtal devices.

Dixtal

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
All Dixtal devices.


Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5


Cable length: 1,0m >




Dixtal

One piece EKG
10 lead cable.

Options		Code
Plug / Banana		BIOECG.DI.P10F



< Cable length: 3,0m >




Compatible with: All Dixtal devices.

Dixtal




< Cable length: 3,0m >




Temperature probe.

Compatible with:
All Dixtal devices.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.400PT	
Skin TITANIUM® 5,5mm		BIOTEMP.400PTP	
Timpanic		BIOTEMP.400T	
Adult esophageal / rectal		BIOTEMP.400RA	
Pediatric esophageal / rectal		BIOTEMP.400RP	



< Cable length: 2,7m >




Dräger / Siemens

SpO2 sensor.

Compatible with:
SC6002XL, 7000, SC8000, SC9000X, Infinity Vista and Infinity Delta.

Options		Code	
Adult Clip		BIOOX.DRS.AC3	
Adult Silicone		BIOOX.DRS.AS3	
Pediatric Clip		BIOOX.DRS.PC3	
Pediatric Silicone		BIOOX.DRS.PS3	
Neonate / Multisite		BIOOX.DRS.Y3	
Wrap silicone		BIOOX.DRS.WS3	



< Cable length: 0,9m >




Dräger / Siemens


SpO2 short sensor.
this product requires an extension cable

Compatible with:
SC6002XL, 7000, SC8000, SC9000X, Infinity Vista and Infinity Delta.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



< Cable length: 2,0m >



Dräger / Siemens


SpO2 intermediary cable.

Compatible with: Infinity Vista, Delta, Gamma, Kappa MultiMed Pod; Siemens SC 6002XL, SC7000, SC8000, SC9000XL.

Description	Code	
SpO2 intermediary cable	BIOOX.DRS.IC2	



<Cable length: 2,7m>










Dräger / Siemens

SpO2 sensor.

Compatible with:

Sirecust 630, Sirecust 700, Sirecust 732, Sirem SpO2 Module, SC 6002 XL, 7000, 8000, 9000 XL.

Options		Code	
Adult Clip		BIOOX.SI.AC3	
Adult Silicone		BIOOX.SI.AS3	
Pediatric Clip		BIOOX.SI.PC3	
Pediatric Silicone		BIOOX.SI.PS3	
Neonate / Multisite		BIOOX.SI.Y3	
Wrap silicone		BIOOX.SI.WS3	



< Cable length: 0,9m >









Dräger / Siemens


SpO2 short sensor.

this product requires an extension cable

Compatible with:

Sirecust 630, Sirecust 700, Sirecust 732, Sirem SpO2 Module, SC 6002 XL, 7000, 8000, 9000 XL.

Options		Code	
Adult Clip		BIOOX.BC.AC1	
Adult Silicone		BIOOX.BC.AS1	
Pediatric Clip		BIOOX.BC.PC1	
Pediatric Silicone		BIOOX.BC.PS1	
Neonate / Multisite		BIOOX.BC.Y1	
Wrap silicone		BIOOX.BC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Dräger / Siemens





SpO2 intermediary cable.


Description	Code	
SpO2 intermediary cable	BIOOX.SI.IC2	

Compatible with: Sirecust 630, Sirecust 700, Sirecust 732, Sirem SpO2 Module, SC 6002 XL, 7000, 8000, 9000 XL.


Dräger / Siemens

One piece ECG cable
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.DR.S3F	
3 leads clamp		BIOECG.DR.C3F	
5 leads snap		BIOECG.DR.S5F	
5 leads clamp		BIOECG.DR.C5F	



< Cable length: 3,5m >







Compatible with:

UM1, UM3, UM3.1, RM1, PM8010, PM8014, Babyguard 8000, Cato, PB8800, PM8060; Series 8000, Cicero, Parameterbox.


Dräger / Siemens

ECG trunk cable with leads
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.DR.S3D	
3 leads clamp		BIOECG.DR.C3D	
5 leads snap		BIOECG.DR.S5D	
5 leads clamp		BIOECG.DR.C5D	



< Cable length: 4,0m >






Compatible with:

UM1, UM3, UM3.1, RM1, PM8010, PM8014, Babyguard 8000, Cato, PB8800, PM8060; Series 8000, Cicero, Parameterbox.


Dräger / Siemens

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.DR.3T	
5 leads		BIOECG.DR.5T	



< Cable length: 3,0m >








Compatible with: UM1, UM3, UM3.1, RM1, PM8010, PM8014, Babyguard 8000, Cato, PB8800, PM8060; Series 8000, Cicero, Parameterbox.

Dräger / Siemens


One piece ECG cable
(3 or 5 leads)

Compatible with:
Vitalguard 2000.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	







< Cable length: 3,5m >




Dräger / Siemens


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
Vitalguard 2000.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




Dräger / Siemens

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	



< Cable length: 3,0m >







Compatible with: Vitalguard 2000.

Dräger / Siemens

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
All Dräger / Siemens devices.

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	



Cable length: 1,0m >



Dräger / Siemens


Temperature probe.

Compatible with: Serie 8000.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.2001PT	
Skin TITANIUM® 5,5mm		BIOTEMP.2001PTP	





< Cable length: 1,5m >



Dräger / Siemens

Temperature probe.

Compatible with: Dräger Caleo Series.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.2200PT	
Skin TITANIUM® 5,5mm		BIOTEMP.2200PTP	








< Cable length: 1,5m >



Dräger / Siemens


Temperature probe.

Compatible with:
Infinity Vista, Delta, Gamma, Kappa.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.1100PT	
Skin TITANIUM® 5,5mm		BIOTEMP.1100PTP	
Timpanic		BIOTEMP.1100T	
Adult esophageal / rectal		BIOTEMP.1100RA	
Pediatric esophageal / rectal		BIOTEMP.1100RP	




< Cable length: 3,0m >





< Cable length: 2,7m >



Ecafix


SpO2 sensor.

Compatible with:
Active ES, Active ESP.

Options		Code	
Adult Clip		BIOOX.PA.AC3	
Adult Silicone		BIOOX.PA.AS3	
Pediatric Clip		BIOOX.PA.PC3	
Pediatric Silicone		BIOOX.PA.PS3	
Neonate / Multisite		BIOOX.PA.Y3	
Wrap silicone		BIOOX.PA.WS3	



< Cable length: 2,7m >



Ecafix

SpO2 sensor.


Compatible with:
Active S, Pluma M, Pluma P. Oxi Cap.

Options		Code	
Adult Clip		BIOOX.EC.AC3	
Adult Silicone		BIOOX.EC.AS3	
Pediatric Clip		BIOOX.EC.PC3	
Pediatric Silicone		BIOOX.EC.PS3	
Neonate / Multisite		BIOOX.EC.Y3	
Wrap silicone		BIOOX.EC.WS3	

Ecafix

One piece ECG cable
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.EC.S3F	
3 leads clamp		BIOECG.EC.C3F	
5 leads snap		BIOECG.EC.S5F	
5 leads clamp		BIOECG.EC.C5F	



< Cable length: 3,5m >








Compatible with: All Ecafix devices.

Ecafix

ECG trunk cable with leads
(3 or 5 leads)



Compatible with:
All Ecafix devices.


Options		Code	
3 leads snap		BIOECG.EC.S3D	
3 leads clamp		BIOECG.EC.C3D	
5 leads snap		BIOECG.EC.S5D	
5 leads clamp		BIOECG.EC.C5D	



Ecafix

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.EC.3T	
5 leads		BIOECG.EC.5T	







Compatible with: All Ecafix devices.


Ecafix

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
All Ecafix devices.


Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >





< Cable length: 2,7m >










Edan

SpO2 sensor.

Compatible with:

M3, M3B, M8, M8B, M9, M99A.

Options		Code	
Adult Clip		BIOOX.ED.AC3	
Adult Silicone		BIOOX.ED.AS3	
Pediatric Clip		BIOOX.ED.PC3	
Pediatric Silicone		BIOOX.ED.PS3	
Neonate / Multisite		BIOOX.ED.Y3	
Wrap silicone		BIOOX.ED.WS3	



< Cable length: 0,9m >









Edan


SpO2 short sensor.

this product requires an extension cable

Compatible with:

M3, M3B, M8, M8B, M9, M99A.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Edan


SpO2 intermediary cable.


Compatible with: M3, M3B, M8, M8B, M9, M99A.

Description	Code	
SpO2 intermediary cable	BIOOX.ED.IC2	


Edan

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.ED.P10F	




< Cable length: 3,0m >




Compatible with: All Edan devices.


Electronica Trentina / Cardiette

Options		Code	
Plug / Banana		BIOECG.NK.P10F	

One piece EKG 10 lead cable.



< Cable length: 3,0m >




Compatible with:

Autorouler 12/1, Autorouler 12/3, Microruler 12/1, Microruler 12/3, Microtel, Excel 101, Excel 103, Excel 106.



< Cable length: 2,7m >










Emai / Transmai

SpO2 sensor.

Compatible with:

OXP-10, MX 100, MX 300, PX 200, 3303, Capnocheck 8400, 9004.

Options		Code	
Adult Clip		BIOOX.EM.AC3	
Adult Silicone		BIOOX.EM.AS3	
Pediatric Clip		BIOOX.EM.PC3	
Pediatric Silicone		BIOOX.EM.PS3	
Neonate / Multisite		BIOOX.EM.Y3	
Wrap silicone		BIOOX.EM.WS3	



< Cable length: 0,9m >









Emai / Transmai


SpO2 short sensor.

this product requires an extension cable

Compatible with:

RX100, RX300.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection



Emai / Transmai

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.EM.IC2	


Compatible with: OXP-10, MX 100, MX 300, PX 200, 3303, Capnocheck 8400, 9004.

Emai / Transmai

One piece ECG cable
(3 or 5 leads)

Compatible with:
RX 300 y MX series.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



< Cable length: 3,5m >




Emai / Transmai


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
RX 300 y MX series.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	




< Cable length: 4,0m >



Emai / Transmai

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	



< Cable length: 3,0m >







Compatible with: RX 300 y MX series.

Emai / Transmai

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
RX 300 y MX series.

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	



Cable length: 1,0m >




Esaote

One piece ECG
cable 3 leads only.

Compatible with: All Esaote devices.

Options		Code	
3 leads snap		BIOECG.ES.S3F	
3 leads clamp		BIOECG.ES.C3F	





< Cable length: 3,5m >




Esaote

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.SC.P10F	




< Cable length: 3,0m >



Compatible with: EKG P120.


Esaote

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.SC2.P10F	



< Cable length: 3,0m >



Compatible with: EKG P80, P8000 Series.

Fanem

Temperature probe.



Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.1000PT	
Skin TITANIUM® 5,5mm		BIOTEMP.1000PTP	

< Cable length: 1,5m >




Compatible with: All Fanem devices.

Feas


One piece ECG cable
(3 or 5 leads)

Compatible with:
All Feas devices.

Options		Code	
3 leads snap		BIOECG.FE.S3F	
3 leads clamp		BIOECG.FE.C3F	
5 leads snap		BIOECG.FE.S5F	
5 leads clamp		BIOECG.FE.C5F	



< Cable length: 3,5m >




Feas


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
All Feas devices.

Options		Code	
3 leads snap		BIOECG.FE.S3D	
3 leads clamp		BIOECG.FE.C3D	
5 leads snap		BIOECG.FE.S5D	
5 leads clamp		BIOECG.FE.C5D	




< Cable length: 4,0m >




Feas

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.FE.3T	
5 leads		BIOECG.FE.5T	



< Cable length: 3,0m >







Compatible with: All Feas devices.

Feas

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
All Feas devices.


Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >



Fischer & Paykel

Temperature probe.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.900PT	




< Cable length: 1,5m >



Compatible with: IW 700, IW 800 & IW 900 Series.

Fischer & Paykel

Temperature probe Dual.

Options		Code	
Air		BIOTEMP.900AD	



< Cable length: 2,0m >



Compatible with: MR700.









< Cable length: 2,7m >




Fukuda Denshi

SpO2 sensor.

Compatible with:
HG302 Module, DS-5100E, DS-5300W.

Options		Code	
Adult Clip		BIOOX.NC.AC3	
Adult Silicone		BIOOX.NC.AS3	
Pediatric Clip		BIOOX.NC.PC3	
Pediatric Silicone		BIOOX.NC.PS3	
Neonate / Multisite		BIOOX.NC.Y3	
Wrap silicone		BIOOX.NC.WS3	



< Cable length: 0,9m >




Fukuda Denshi

SpO2 short sensor.
this product requires an extension cable

Compatible with:
HG302 Module, DS-5100E, DS-5300W.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Fukuda Denshi

SpO2 intermediary cable.


Compatible with: HG302 Module, DS-5100E, DS-5300W.

Description	Code	
SpO2 intermediary cable	BIOOX.NC.IC2	


Fukuda Denshi

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.FU.P10F	



< Cable length: 3,0m >




Compatible with: FX-101, FX-102F, FX-3010, FX-328UY, FCP 2155, FCP 2201, FCP 4101, FCP 4101A, CardiMax FX-7101, FX-7202, FX-7402.

GE - Critikon


One piece ECG cable
(3 or 5 leads)

Compatible with:
Critikon 8100T, Dinamap, 8700,
8710, 8720, 9700, 9710, 9720.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



< Cable length: 3,5m >








GE - Critikon

ECG trunk cable with leads
(3 or 5 leads)


Compatible with:

Critikon 8100T, Dinamap, 8700,
8710, 8720, 9700, 9710, 9720.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




GE - Critikon

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	



< Cable length: 3,0m >







Compatible with: Critikon 8100T, Dinamap, 8700, 8710, 8720, 9700, 9710, 9720.

GE - Critikon

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:

Critikon 8100T, Dinamap, 8700, 8710,
8720, 9700, 9710, 9720.

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >





< Cable length: 2,7m >




GE - Datex Ohmeda

SpO2 sensor.


Compatible with:

AS/3, AS/5, CS/3 Modular Monitors with M-series modules, CH/S, CH/2S, CH/RS, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Cardiocap S/5 Light Monitor, Capnomac Ultima, M-Estpr, Oscar II, Satellite, Satellite Plus II, Trans Plus, ULTSVL, ULTSVO, ULTSIO, ULTS00, OSP-200.

Options		Code	
Adult Clip		BIOOX.DA.AC3	
Adult Silicone		BIOOX.DA.AS3	
Pediatric Clip		BIOOX.DA.PC3	
Pediatric Silicone		BIOOX.DA.PS3	
Neonate / Multisite		BIOOX.DA.Y3	
Wrap silicone		BIOOX.DA.WS3	



< Cable length: 2,7m >




GE - Datex Ohmeda

SpO2 sensor.


Compatible with:

S/5, i4 Modular Monitors con E-Series, Corometrics 259 with Trusignal Ohmeda® Technology.

Options		Code	
Adult Clip		BIOOX.MQO.AC3	
Adult Silicone		BIOOX.MQO.AS3	
Pediatric Clip		BIOOX.MQO.PC3	
Pediatric Silicone		BIOOX.MQO.PS3	
Neonate / Multisite		BIOOX.MQO.Y3	
Wrap silicone		BIOOX.MQO.WS3	



< Cable length: 2,7m >



GE - Datex Ohmeda

SpO2 sensor.


Compatible with:

Ohmeda Trusat.

Options		Code	
Adult Clip		BIOOX.GET.AC3	
Adult Silicone		BIOOX.GET.AS3	
Pediatric Clip		BIOOX.GET.PC3	
Pediatric Silicone		BIOOX.GET.PS3	
Neonate / Multisite		BIOOX.GET.Y3	
Wrap silicone		BIOOX.GET.WS3	



< Cable length: 2,7m >



GE - Datex Ohmeda

SpO2 sensor.

Compatible with:

Ohmeda 3700, 3710, 3770, 3775, 3800, 3900, Biox 3740, 3760, 3800, 4700 Oxicap, 5250 RGM, Sac II, Oxi-Tip Plus, Tuffsat, S/5.

Options		Code	
Adult Clip		BIOOX.OHT.AC3	
Adult Silicone		BIOOX.OHT.AS3	
Pediatric Clip		BIOOX.OHT.PC3	
Pediatric Silicone		BIOOX.OHT.PS3	
Neonate / Multisite		BIOOX.OHT.Y3	
Wrap silicone		BIOOX.OHT.WS3	



< Cable length: 0,9m >



GE - Datex Ohmeda

SpO2 short sensor.

this product requires an extension cable

Compatible with:

Datex Ohmeda Trusignal devices.

Options		Code	
Adult Clip		BIOOX.GES.AC1	
Adult Silicone		BIOOX.GES.AS1	
Pediatric Clip		BIOOX.GES.PC1	
Pediatric Silicone		BIOOX.GES.PS1	
Neonate / Multisite		BIOOX.GES.Y1	
Wrap silicone		BIOOX.GES.WS1	



< Cable length: 0,9m >



GE - Datex Ohmeda

SpO2 short sensor.

this product requires an extension cable

Compatible with:

Universal Oxytip.

Options		Code	
Adult Clip		BIOOX.GEO.AC1	
Adult Silicone		BIOOX.GEO.AS1	
Pediatric Clip		BIOOX.GEO.PC1	
Pediatric Silicone		BIOOX.GEO.PS1	
Neonate / Multisite		BIOOX.GEO.Y1	
Wrap silicone		BIOOX.GEO.WS1	



< Cable length: 0,9m >




GE - Datex Ohmeda

SpO2 short sensor.
this product requires an extension cable

Compatible with:
AS/3, CS/3, CH/S, CH/2S, CH-RS, SAF-F, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Capnomac Ultima, Light Monitor, M-Estpr, Oscar II, OSP-200, Satellite Trans, Satellite, Satellite Plus II, Trans Plus, System 5, ULTS00, ULTSIO, ULTSVI, ULTSVO.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

GE - Datex Ohmeda

SpO2 intermediary cable.


Compatible with: AS/3, AS/5, CS/3 Modular Monitors with M-series modules, CH/S, CH/2S, CH/RS, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Cardiocap S/5 Light Monitor, Capnomac Ultima, M-Estpr, Oscar II, Satellite, Satellite Plus II, Trans Plus, ULTSVL, ULTSVO, ULTSIO, ULTS00, OSP-200.

Description	Code	
SpO2 intermediary cable	BIOOX.DA.IC2	



sensor connection

< Cable length: 2,0m >



device connection

GE - Datex Ohmeda

SpO2 intermediary cable.

Compatible with: Ohmeda 3700, 3710, 3770, 3775, 3800, 3900, Biox 3740, 3760, 3800, 4700 Oxicap, 5250 RGM, Sac II, Oxi-Tip Plus, Tuffsat, S/5.

Description	Code	
SpO2 intermediary cable	BIOOX.GEO.IC2	99.80



< Cable length: 2,0m >



device connection

GE - Datex Ohmeda

SpO2 intermediary cable.


Compatible with: CARESCAPE™ modular monitors and S/5, AS/3 and CS/3 modular monitors with M-NE12ST(P)R, M-(N)EST(P)(R) or M-(P)RESTN modules; Cardiocap/5 and S/5 Light Monitors.

Description	Code	
SpO2 intermediary cable	BIOOX.DT.IC2	



sensor connection

< Cable length: 2,0m >



device connection

GE - Datex Ohmeda

SpO2 intermediary cable.

Compatible with: Ohmeda Trusat.

Description	Code	
SpO2 intermediary cable	BIOOX.GET.IC2	



sensor connection

< Cable length: 2,0m >



device connection

GE - Datex Ohmeda

SpO2 intermediary cable.


Compatible with: CARESCAPE™ modular monitors and S/5, AS/3 and CS/3 modular monitors with E-(P)RESTN or E-PSM(P) modules. CARESCAPE™ V100 with TruSignal® option, Responder 2000.

Description	Code	
SpO2 intermediary cable	BIOOX.MQO.IC2	



sensor connection

< Cable length: 2,0m >



device connection

GE - Datex Ohmeda





SpO2 intermediary cable.


Compatible with: Handheld pulse oximeters (GE Ohmeda TuffSat®, 3775), tabletop pulse oximeters (3700 series, 3800, 3900, 3900P), respiratory gas monitors (4700, 5250 RGM), CARESCAPE modular monitors and S/5, AS/3, CS/3 modular monitors with M-OSAT module, Cardiocap™/5 with N-XOSAT option.

Description	Code	
SpO2 intermediary cable	BIOOX.OHT.IC2	


GE - Datex Ohmeda

One piece ECG cable (3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.DA.S3F	
3 leads clamp		BIOECG.DA.C3F	
5 leads snap		BIOECG.DA.S5F	
5 leads clamp		BIOECG.DA.C5F	







< Cable length: 3,5m >




Compatible with: AS/3, CS/3, S/5, Cardiocap I, Cardiocap II, Cardiocap 5, Light Monitor, Satlite plus.

GE - Datex Ohmeda

ECG trunk cable with leads
(3 or 5 leads)



Options		Code	
3 leads snap		BIOECG.DA.S3D	
3 leads clamp		BIOECG.DA.C3D	
5 leads snap		BIOECG.DA.S5D	
5 leads clamp		BIOECG.DA.C5D	




Compatible with: AS/3, CS/3, S/5, Cardiocap I, Cardiocap II, Cardiocap 5, Light Monitor, Satlite plus.

GE - Datex Ohmeda

ECG trunk cable
(3 or 5 leads).
this product requires a lead set





Options		Code	
3 leads		BIOECG.DA.3T	
5 leads		BIOECG.DA.5T	



Compatible with: AS/3, CS/3, S/5, Cardiocap I, Cardiocap II, Cardiocap 5, Light Monitor, Satlite plus.

GE - Datex Ohmeda

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	



Compatible with:
AS/3, CS/3, S/5, Cardiocap I, Cardiocap II,
Cardiocap 5, Light Monitor, Satlite plus.


Cable length: 1,0m >



GE - Datex Ohmeda

Temperature probe.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.300PT	
Skin TITANIUM® 5,5mm		BIOTEMP.300PTP	




< Cable length: 1,5m >




Compatible with: Care Plus, GC, IC, 3000, 4000, 5000 Series Infant Warmer Systems.


GE - Datex Ohmeda

Temperature probe.

Options		Code	
Air		BIOTEMP.300A3	



< Cable length: 1,5m >



Compatible with: All Ohmeda Care plus devices.


GE - Datex Ohmeda

Temperature probe.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.1700PT	



< Cable length: 1,5m >




Compatible with: All Ohmeda Giraffe devices.

GE - Datex Ohmeda






Temperature probe.



< Cable length: 3,0m >




Compatible with:
All GE - Datex Ohmeda devices.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.400PT	
Skin TITANIUM® 5,5mm		BIOTEMP.400PTP	
Timpanic		BIOTEMP.400T	
Adult esophageal / rectal		BIOTEMP.400RA	
Pediatric esophageal / rectal		BIOTEMP.400RP	



< Cable length: 2,7m >










GE - Marquette

SpO2 sensor.


Compatible with:

EAGLE, DASH, SOLAR Series, GE Responder 2000, Responder 3000, TRAM Modular x50SL Series Nellcor® Technology.

Options		Code	
Adult Clip		BIOOX.MQ.AC3	
Adult Silicone		BIOOX.MQ.AS3	
Pediatric Clip		BIOOX.MQ.PC3	
Pediatric Silicone		BIOOX.MQ.PS3	
Neonate / Multisite		BIOOX.MQ.Y3	
Wrap silicone		BIOOX.MQ.WS3	



< Cable length: 2,7m >










GE - Marquette

SpO2 sensor.

Compatible with:

Dash 2000, 3000, 4000, 5000, i/9500 (Tram 451N5/851N5). Solar 7000, 8000 e Eagle 1000, 3000, Smk SMV, Corometrics 250 with Masimo technology.

Options		Code	
Adult Clip		BIOOX.MQM.AC3	
Adult Silicone		BIOOX.MQM.AS3	
Pediatric Clip		BIOOX.MQM.PC3	
Pediatric Silicone		BIOOX.MQM.PS3	
Neonate / Multisite		BIOOX.MQM.Y3	
Wrap silicone		BIOOX.MQM.WS3	










< Cable length: 2,7m >

GE - Marquette

SpO2 sensor.

Compatible with:
Corometrics 250cx and Marquette monitors with Oximax® technology.

Options		Code	
Adult Clip		BIOOX.MQX.AC3	
Adult Silicone		BIOOX.MQX.AS3	
Pediatric Clip		BIOOX.MQX.PC3	
Pediatric Silicone		BIOOX.MQX.PS3	
Neonate / Multisite		BIOOX.MQX.Y3	
Wrap silicone		BIOOX.MQX.WS3	









< Cable length: 2,7m >

GE - Marquette

SpO2 sensor.

Compatible with:
Hellige Eagle 1000, 3000, SMU610, 611, SMK210, SMV, Series with Nellcor® technology (standard).

Options		Code	
Adult Clip		BIOOX.GEH.AC3	
Adult Silicone		BIOOX.GEH.AS3	
Pediatric Clip		BIOOX.GEH.PC3	
Pediatric Silicone		BIOOX.GEH.PS3	
Neonate / Multisite		BIOOX.GEH.Y3	
Wrap silicone		BIOOX.GEH.WS3	









< Cable length: 0,9m >


GE - Marquette

SpO2 short sensor.

this product requires an extension cable


Compatible with:
All GE - Marquette devices with Nellcor® technology.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >




device connection

GE - Marquette

SpO2 intermediary cable.


Description	Code	
SpO2 intermediary cable	BIOOX.MQ.IC2	

Compatible with: EAGLE, DASH, SOLAR Series, GE Responder 2000, Responder 3000, TRAM Modular x50SL Series with Nellcor® Technology.



sensor connection

< Cable length: 2,0m >




device connection

GE - Marquette

SpO2 intermediary cable.


Description	Code	
SpO2 intermediary cable	BIOOX.MQM.IC2	

Compatible with: Dash 2000, 3000, 4000, 5000, i/9500 (Tram 451N5/851N5). Solar 7000, 8000 e Eagle 1000, 3000, Smk SMV, Corometrics 250 with Masimo® technology.



sensor connection

< Cable length: 2,0m >



device connection

GE - Marquette





SpO2 intermediary cable.


Description	Code	
SpO2 intermediary cable	BIOOX.GEH.IC2	

Compatible with: Hellige Eagle 1000, 3000, SMU610, 611, SMK210, SMV, Series with Nellcor® technology (standard).

GE - Marquette

One piece ECG cable
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.MQ.S3F	
3 leads clamp		BIOECG.MQ.C3F	
5 leads snap		BIOECG.MQ.S5F	
5 leads clamp		BIOECG.MQ.C5F	







< Cable length: 3,5m >




Compatible with: Dash 1000, Dash Pro 2000, Dash 2500, Dash Pro 3000, Dash Pro 4000, Dash 5000, Eagle Systems, MAC-Lab, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems and Carescape.

GE - Marquette

ECG trunk cable with leads
(3 or 5 leads)



Options		Code	
3 leads snap		BIOECG.MQ.S3D	
3 leads clamp		BIOECG.MQ.C3D	
5 leads snap		BIOECG.MQ.S5D	
5 leads clamp		BIOECG.MQ.C5D	




Compatible with: Dash 1000, Dash Pro 2000, Dash 2500, Dash Pro 3000, Dash Pro 4000, Dash 5000, Eagle Systems, MAC-Lab, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems and Careescape.

GE - Marquette

ECG trunk cable
(3 or 5 leads).
this product requires a lead set





Options		Code	
3 leads		BIOECG.MQ.3T	
5 leads		BIOECG.MQ.5T	



Compatible with: Dash 1000, Dash Pro 2000, Dash 2500, Dash Pro 3000, Dash Pro 4000, Dash 5000, Eagle Systems, MAC-Lab, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems and Careescape.

GE - Marquette

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >







Compatible with: Dash 1000, Dash Pro 2000, Dash 2500, Dash Pro 3000, Dash Pro 4000, Dash 5000, Eagle Systems, MAC-Lab, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems and Careescape.


GE - Marquette

One piece ECG cable (3 or 5 leads)


Compatible with:

Hellige, Servomed 4. Generation, SCP 851, SCP 852, SCP 844, SCS 804, Cardioserv, VICOM-SM.

Options		Code	
3 leads snap		BIOECG.MQH.S3F	
3 leads clamp		BIOECG.MQH.C3F	
5 leads snap		BIOECG.MQH.S5F	
5 leads clamp		BIOECG.MQH.C5F	



< Cable length: 3,5m >








GE - Marquette

ECG trunk cable with leads (3 or 5 leads)

Compatible with:

Hellige, Servomed 4. Generation, SCP 851, SCP 852, SCP 844, SCS 804, Cardioserv, VICOM-SM.

Options		Code	
3 leads snap		BIOECG.MQH.S3D	
3 leads clamp		BIOECG.MQH.C3D	
5 leads snap		BIOECG.MQH.S5D	
5 leads clamp		BIOECG.MQH.C5D	





< Cable length: 4,0m >




GE - Marquette

ECG trunk cable (3 or 5 leads).

this product requires a lead set

Options		Code	
3 leads		BIOECG.MQH.3T	
5 leads		BIOECG.MQH.5T	



< Cable length: 3,0m >



Compatible with: Hellige, Servomed 4. Generation, SCP 851, SCP 852, SCP 844, SCS 804, Cardioserv, VICOM-SM.





GE - Marquette

ECG lead set (3 or 5 leads)

this product requires a trunk cable

Compatible with:

Hellige, Servomed 4. Generation, SCP 851, SCP 852, SCP 844, SCS 804, Cardioserv, VICOM-SM.


Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	


Cable length: 1,0m >




GE - Marquette

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.HE.P10F	




< Cable length: 3,0m >




Compatible with: EK-36, EK-41, EK-43, EK-53, EK-56, EK-403, EK-413, EK-512.


GE - Marquette

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.HE.P10F	



< Cable length: 3,0m >




Compatible with: MAC 400, MAC 500, MAC 1100, MAC 1200.






GE - Marquette

Temperature
probe.




< Cable length: 3,0m >



Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.1200PT	
Skin TITANIUM® 5,5mm		BIOTEMP.1200PTP	
Timpanic		BIOTEMP.1200T	
Adult esophageal / rectal		BIOTEMP.1200RA	
Pediatric esophageal / rectal		BIOTEMP.1200RP	

Compatible with:

Carescape, Dash, Eagle,
Tram, Solar, Mac- Lab,
Cath-Lab.









< Cable length: 2,7m >




General Meditech

SpO2 sensor.

Options		Code	
Adult Clip		BIOOX.MT.AC3	
Adult Silicone		BIOOX.MT.AS3	
Pediatric Clip		BIOOX.MT.PC3	
Pediatric Silicone		BIOOX.MT.PS3	
Neonate / Multisite		BIOOX.MT.Y3	
Wrap silicone		BIOOX.MT.WS3	

Compatible with:

G2A, G3A, G3C, G3D,
G3F, G3G, G3H, G3L.









< Cable length: 2,7m >



General Meditech

SpO2 sensor.

Compatible with:
G1B.


Options		Code
Adult Clip		BIOOX.MTD.AC3
Adult Silicone		BIOOX.MTD.AS3
Pediatric Clip		BIOOX.MTD.PC3
Pediatric Silicone		BIOOX.MTD.PS3
Neonate / Multisite		BIOOX.MTD.Y3
Wrap silicone		BIOOX.MTD.WS3

Innomed


One piece EKG
10 lead cable.



< Cable length: 3,0m >



Compatible with: Heart Mirror 1 IKO - 3 - 3D, Heart screen 80G - 112C - 112D.

Options		Code
Plug / Banana		BIOECG.SC.P10F









< Cable length: 2,7m >



Instramed

SpO2 sensor.





Compatible with:
MinisCodepe II, Cardioversor
HS 03 (DB-9 s/ext).


Options		Code
Adult Clip		BIOOX.BC.AC3
Adult Silicone		BIOOX.BC.AS3
Pediatric Clip		BIOOX.BC.PC3
Pediatric Silicone		BIOOX.BC.PS3
Neonate / Multisite		BIOOX.BC.Y3
Wrap silicone		BIOOX.BC.WS3

Instramed


One piece ECG cable
(3 or 5 leads)

Compatible with:
Inmax Series.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	







< Cable length: 3,5m >




Instramed


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
Inmax Series.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




Instramed

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	



< Cable length: 3,0m >







Compatible with: Inmax Series.

Instramed


ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
Inmax Series.

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >





< Cable length: 2,7m >



Invivo Research

SpO2 sensor.


Compatible with:

4500, 4500 Plus, Millenia, Sat Scout Omega 1445, Omni-Trak TVS.

Options		Code	
Adult Clip		BIOOX.IR.AC3	
Adult Silicone		BIOOX.IR.AS3	
Pediatric Clip		BIOOX.IR.PC3	
Pediatric Silicone		BIOOX.IR.PS3	
Neonate / Multisite		BIOOX.IR.Y3	
Wrap silicone		BIOOX.IR.WS3	



< Cable length: 0,9m >



Invivo Research

SpO2 short sensor.

this product requires an extension cable

Compatible with:

All Invivo Research devices with DB9 6 pin connector.

Options		Code	
Adult Clip		BIOOX.IRD.AC1	
Adult Silicone		BIOOX.IRD.AS1	
Pediatric Clip		BIOOX.IRD.PC1	
Pediatric Silicone		BIOOX.IRD.PS1	
Neonate / Multisite		BIOOX.IRD.Y1	
Wrap silicone		BIOOX.IRD.WS1	



< Cable length: 0,9m >



Invivo Research


SpO2 short sensor.

this product requires an extension cable

Compatible with:

4500 Plus 2 with Nellcor® technology.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



< Cable length: 2,0m >



Invivo Research

sensor connection

device connection


SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.IR.IC2	

Compatible with: 4500, 4500 Plus, Millenia, Sat Scout Omega 1445, Omni-Trak TVS.









< Cable length: 2,7m >



Ivy Biomedical


SpO2 sensor.


Compatible with:
All Ivy Biomedical devices.

Options		Code	
Adult Clip		BIOOX.OHT.AC3	
Adult Silicone		BIOOX.OHT.AS3	
Pediatric Clip		BIOOX.OHT.PC3	
Pediatric Silicone		BIOOX.OHT.PS3	
Neonate / Multisite		BIOOX.OHT.Y3	
Wrap silicone		BIOOX.OHT.WS3	


Kenz

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.CA.P10F	



< Cable length: 3,0m >



Compatible with: 103, 106, 107, 302, 303, 1201.









< Cable length: 2,7m >



Kontron

SpO2 sensor.

Compatible with:
7840, 7845, Minimon 7137, Micromon,
Supermon, Kolormon.

Options		Code	
Adult Clip		BIOOX.KO.AC3	
Adult Silicone		BIOOX.KO.AS3	
Pediatric Clip		BIOOX.KO.PC3	
Pediatric Silicone		BIOOX.KO.PS3	
Neonate / Multisite		BIOOX.KO.Y3	
Wrap silicone		BIOOX.KO.WS3	



< Cable length: 0,9m >




Kontron

SpO2 short sensor.
this product requires an extension cable

Compatible with:
7840, 7845, Minimon 7137, Micromon,
Supermon, Kolormon.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Kontron

SpO2 intermediary cable.

Compatible with: 7840, 7845, Minimon 7137, Micromon, Supermon, Kolormon.

Description	Code	
SpO2 intermediary cable	BIOOX.KO.IC2	

Kontron



One piece EKG
10 lead cable.




Compatible with: EKG P80.

Options		Code	
Plug / Banana		BIOECG.SC.P10F	

< Cable length: 3,0m >



< Cable length: 2,7m >




Ktmed Epox

SpO2 sensor.

Compatible with:
KTPE-01.

Options		Code	
Adult Clip		BIOOX.EP.AC3	
Adult Silicone		BIOOX.EP.AS3	
Pediatric Clip		BIOOX.EP.PC3	
Pediatric Silicone		BIOOX.EP.PS3	
Neonate / Multisite		BIOOX.EP.Y3	
Wrap silicone		BIOOX.EP.WS3	



< Cable length: 0,9m >




Ktmed Epox

SpO2 short sensor.
this product requires an extension cable

Compatible with:
KTPE-01.


Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	

Ktmed Epox



sensor connection

< Cable length: 2,0m >



device connection

SpO2 intermediary cable.

Compatible with: KTPE-01.

Description	Code	
SpO2 intermediary cable	BIOOX.EP.IC2	



< Cable length: 2,7m >




MDE

SpO2 sensor.

Compatible with:
Escort 1 / E 300, EsCodert 100 OPT
11A, EsCodert 300A, Escort II OPT30,
Prism, 20300, 20401, 20403.

Options		Code	
Adult Clip		BIOOX.NC.AC3	
Adult Silicone		BIOOX.NC.AS3	
Pediatric Clip		BIOOX.NC.PC3	
Pediatric Silicone		BIOOX.NC.PS3	
Neonate / Multisite		BIOOX.NC.Y3	
Wrap silicone		BIOOX.NC.WS3	



< Cable length: 2,7m >



Medexel


SpO2 sensor.

Compatible with:
All Medexel devices with DB9
connector.

Options		Code	
Adult Clip		BIOOX.NC.AC3	
Adult Silicone		BIOOX.NC.AS3	
Pediatric Clip		BIOOX.NC.PC3	
Pediatric Silicone		BIOOX.NC.PS3	
Neonate / Multisite		BIOOX.NC.Y3	
Wrap silicone		BIOOX.NC.WS3	



Medix

Temperature probe.


Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.800PT	
Skin TITANIUM® 5,5mm		BIOTEMP.800PTP	



< Cable length: 1,5m >







**Compatible with:** All Medix devices.

< Cable length: 2,7m >

**Medinova**


SpO2 sensor.

Compatible with:
All Medinova devices.

Options		Code	
Adult Clip		BIOOX.MV.AC3	
Adult Silicone		BIOOX.MV.AS3	
Pediatric Clip		BIOOX.MV.PC3	
Pediatric Silicone		BIOOX.MV.PS3	
Neonate / Multisite		BIOOX.MV.Y3	
Wrap silicone		BIOOX.MV.WS3	









< Cable length: 2,7m >

**Medlab**


SpO2 sensor.

Compatible with:
Pox 10.

Options		Code	
Adult Clip		BIOOX.MB.AC3	
Adult Silicone		BIOOX.MBAS3	
Pediatric Clip		BIOOX.MB.PC3	
Pediatric Silicone		BIOOX.MB.PS3	
Neonate / Multisite		BIOOX.MB.Y3	
Wrap silicone		BIOOX.MB.WS3	









< Cable length: 2,7m >

**MEK**

SpO2 sensor.

Compatible with:

MP-100, MP-110, MP-400, MP-500, MP-600, MP-1000.

Options		Code	
Adult Clip		BIOOX.MK.AC3	
Adult Silicone		BIOOX.MK.AS3	
Pediatric Clip		BIOOX.MK.PC3	
Pediatric Silicone		BIOOX.MK.PS3	
Neonate / Multisite		BIOOX.MK.Y3	
Wrap silicone		BIOOX.MK.WS3	









< Cable length: 2,7m >


**MEK**

SpO2 sensor.


Compatible with:

MP110 Plus, MP111.

Options		Code	
Adult Clip		BIOOX.NC.AC3	
Adult Silicone		BIOOX.NC.AS3	
Pediatric Clip		BIOOX.NC.PC3	
Pediatric Silicone		BIOOX.NC.PS3	
Neonate / Multisite		BIOOX.NC.Y3	
Wrap silicone		BIOOX.NC.WS3	









< Cable length: 2,7m >

**MEK**

SpO2 sensor.


Compatible with:

PAMO I, PAMO II, PAMO III, MP 570T, MP 700, MP 800C.

Options		Code	
Adult Clip		BIOOX.MKP.AC3	
Adult Silicone		BIOOX.MKP.AS3	
Pediatric Clip		BIOOX.MKP.PC3	
Pediatric Silicone		BIOOX.MKP.PS3	
Neonate / Multisite		BIOOX.MKP.Y3	
Wrap silicone		BIOOX.MKP.WS3	









< Cable length: 2,7m >



MEK


SpO2 sensor.

Compatible with:
MP 1000NT, MP 1300.

Options		Code	
Adult Clip		BIOOX.MK2.AC3	
Adult Silicone		BIOOX.MK2.AS3	
Pediatric Clip		BIOOX.MK2.PC3	
Pediatric Silicone		BIOOX.MK2.PS3	
Neonate / Multisite		BIOOX.MK2.Y3	
Wrap silicone		BIOOX.MK2.WS3	



< Cable length: 0,9m >










MEK

SpO2 short sensor.

this product requires an extension cable

Compatible with:
MP110 Plus, MP111, PAMO I, PAMO II, PAMO III, MP 570T, MP 700, MP 800C, MP 1000NT, MP 1300.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >




device connection

MEK

SpO2 intermediary cable.


Compatible with: MP110 Plus, MP111.

Description	Code	
SpO2 intermediary cable	BIOOX.NC.IC2	



sensor connection

< Cable length: 2,0m >




device connection

MEK

SpO2 intermediary cable.

Compatible with: PAMO I, PAMO II, PAMO III, MP 570T, MP 700, MP 800C.

Description	Code	
SpO2 intermediary cable	BIOOX.MKP.IC2	



sensor connection

< Cable length: 2,0m >



device connection

MEK

SpO2 intermediary cable.





Description	Code	
SpO2 intermediary cable	BIOOX.MK2.IC2	


Compatible with: MP 1000NT, MP 1300.

MEK

One piece ECG cable
(3 or 5 leads)

Compatible with:
1000NT.

Options		Code	
3 leads snap		BIOECG.MK.S3F	
3 leads clamp		BIOECG.MK.C3F	
5 leads snap		BIOECG.MK.S5F	
5 leads clamp		BIOECG.MK.C5F	







< Cable length: 3,5m >



MEK


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
1000NT.

Options		Code	
3 leads snap		BIOECG.MK.S3D	
3 leads clamp		BIOECG.MK.C3D	
5 leads snap		BIOECG.MK.S5D	
5 leads clamp		BIOECG.MK.C5D	






< Cable length: 4,0m >




MEK

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.MK.3T	
5 leads		BIOECG.MK.5T	



< Cable length: 3,0m >







Compatible with: 1000NT.

MEK

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
1000NT.

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >









< Cable length: 2,7m >



Mennen

SpO2 sensor.

Compatible with:
Horizon XL (Nellcor SpO2), Mercury
1200, MR 1300, MR 1303, Envoy.

Options		Code	
Adult Clip		BIOOX.ME.AC3	
Adult Silicone		BIOOX.ME.AS3	
Pediatric Clip		BIOOX.ME.PC3	
Pediatric Silicone		BIOOX.ME.PS3	
Neonate / Multisite		BIOOX.ME.Y3	
Wrap silicone		BIOOX.ME.WS3	









< Cable length: 2,7m >



Mennen

SpO2 sensor.

Compatible with:
Enmove 1000.

Options		Code	
Adult Clip		BIOOX.ME3.AC3	
Adult Silicone		BIOOX.ME3.AS3	
Pediatric Clip		BIOOX.ME3.PC3	
Pediatric Silicone		BIOOX.ME3.PS3	
Neonate / Multisite		BIOOX.ME3.Y3	
Wrap silicone		BIOOX.ME3.WS3	









< Cable length: 0,9m >




Mennen

SpO2 short sensor.
this product requires an extension cable

Compatible with:
Horizon XL (Nellcor SpO2), Mercury
1200, MR 1300, MR 1303, Envoy,
Enmove 1000.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



< Cable length: 2,0m >




Mennen sensor connection

device connection

SpO2 intermediary cable.

Compatible with:
Horizon 2000, Horizon XL (Nellcor SpO2), Mercury 1200, MR 1300, MR 1303, Envoy.

Description	Code	
SpO2 intermediary cable	BIOOX.ME.IC2	



< Cable length: 2,0m >



Mennen sensor connection

device connection

SpO2 intermediary cable.





Compatible with: Enmove 1000.


Description	Code	
SpO2 intermediary cable	BIOOX.ME3.IC2	

Mennen


One piece ECG cable
(3 or 5 leads)

Compatible with:
Horizon 250, 450, 740, 741, 742, 3000, Microsolo, 936 Series, 425/101 Preamp, Cardio-Pak 2000.

Options		Code	
3 leads snap		BIOECG.ME.S3F	
3 leads clamp		BIOECG.ME.C3F	
5 leads snap		BIOECG.ME.S5F	
5 leads clamp		BIOECG.ME.C5F	







< Cable length: 3,5m >




Mennen


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
Horizon 250, 450, 740, 741, 742, 3000, Microsolo, 936 Series, 425/101 Preamp, Cardio-Pak 2000.

Options		Code	
3 leads snap		BIOECG.ME.S3D	
3 leads clamp		BIOECG.ME.C3D	
5 leads snap		BIOECG.ME.S5D	
5 leads clamp		BIOECG.ME.C5D	






< Cable length: 4,0m >




Mennen

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.ME.3T	
5 leads		BIOECG.ME.5T	







< Cable length: 3,0m >



Compatible with: Horizon 250, 450, 740, 741, 742, 3000, Microsolo, 936 Series, 425/101 Preamp, Cardio-Pak 2000.

Mennen

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Compatible with:
Horizon 250, 450, 740, 741, 742,
3000, Microsolo, 936 Series,
425/101 Preamp, Cardio-Pak 2000.





Cable length: 1,0m >




Mennen

One piece ECG
cable (3 or 5 leads)

Compatible with:
Enmove, 1000, 1200,
Envoy, Vitalogik.

Options		Code	
3 leads snap		BIOECG.ME3.S3F	
3 leads clamp		BIOECG.ME3.C3F	
5 leads snap		BIOECG.ME3.S5F	
5 leads clamp		BIOECG.ME3.C5F	



< Cable length: 3,5m >








Mennen

ECG trunk cable with leads
(3 or 5 leads)


Compatible with:

Enmove, 1000, 1200,
Envoy, Vitalogik.

Options		Code
3 leads snap		BIOECG.ME3.S3D
3 leads clamp		BIOECG.ME3.C3D
5 leads snap		BIOECG.ME3.S5D
5 leads clamp		BIOECG.ME3.C5D





< Cable length: 4,0m >




Mennen


ECG trunk cable
(3 or 5 leads).

this product requires a lead set

Options		Code
3 leads		BIOECG.ME3.3T
5 leads		BIOECG.ME3.5T



< Cable length: 3,0m >



Compatible with: Enmove, 1000, 1200, Envoy, Vitalogik.





Mennen

ECG lead set
(3 or 5 leads)

this product requires a trunk cable

Compatible with:

Enmove, 1000, 1200, Envoy,
Vitalogik.

Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5

Cable length: 1,0m >




Mennen






Temperature
probe.


Compatible with:


All Mennen devices with
Plug 1/4 connector.



< Cable length: 3,0m >

Options		Code
Skin TITANIUM® 9,5mm		BIOTEMP.400PT
Skin TITANIUM® 5,5mm		BIOTEMP.400PTP
Timpanic		BIOTEMP.400T
Adult esophageal / rectal		BIOTEMP.400RA
Pediatric esophageal / rectal		BIOTEMP.400RP





< Cable length: 2,7m >



Mindray

SpO2 sensor.


Compatible with:

PM 600, PM 6000 (New), PM 7000, 6201, PM 8000, PM 9000, PM 9000 Express, MEC 100, MEC 200, MEC 1000, MEC 1200, 809B, VS 800; Mennen Enmove 1000, Mercury; Kontron: 7138, 7840, 7845; Medec: MEDEC-1000, M -1000, M- 1100; Creative: PC9000, UP6000; Edan: M3, M3A, M3B, M8, M8A, M8B, M9, M9A, M50, M80, iM8, iM8A, iM8B, iM9, iM9A, iM50, iM80; Goldway, Contec and Danatech, Penlon Intermed PM 8000, Vital Lite, Vital Line, Vital Wave.

Options		Code	
Adult Clip		BIOOX.MI.AC3	
Adult Silicone		BIOOX.MI.AS3	
Pediatric Clip		BIOOX.MI.PC3	
Pediatric Silicone		BIOOX.MI.PS3	
Neonate / Multisite		BIOOX.MI.Y3	
Wrap silicone		BIOOX.MI.WS3	



< Cable length: 2,7m >




Mindray

SpO2 sensor.

Compatible with:

PM 5000, PM 6000, PM 8000.

Options		Code	
Adult Clip		BIOOX.MI2.AC3	
Adult Silicone		BIOOX.MI2.AS3	
Pediatric Clip		BIOOX.MI2.PC3	
Pediatric Silicone		BIOOX.MI2.PS3	
Neonate / Multisite		BIOOX.MI2.Y3	
Wrap silicone		BIOOX.MI2.WS3	



< Cable length: 2,7m >



Mindray

SpO2 sensor.

Compatible with:

BeneView T5, T8 and iPM-8900.

Options		Code	
Adult Clip		BIOOX.MI3.AC3	
Adult Silicone		BIOOX.MI3.AS3	
Pediatric Clip		BIOOX.MI3.PC3	
Pediatric Silicone		BIOOX.MI3.PS3	
Neonate / Multisite		BIOOX.MI3.Y3	
Wrap silicone		BIOOX.MI3.WS3	



< Cable length: 0,9m >



Mindray

SpO2 short sensor.
this product requires an extension cable

Compatible with:
All Mindray devices.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >




device connection

Mindray

SpO2 intermediary cable.


Description	Code	
SpO2 intermediary cable	BIOOX.MI.IC2	

Compatible with: PM 600, PM 6000 (New), PM 7000, 6201, PM 8000, PM 9000, PM 9000 Express, MEC 100, MEC 200, MEC 1000, MEC 1200, 809B, VS 800; Mennen Enmove 1000, Mercury; Kontron: 7138, 7840, 7845; Medec: MEDEC-1000, M-1000, M-1100; Creative: PC9000, UP6000; Edan: M3, M3A, M3B, M8, M8A, M8B, M9, M9A, M50, M80, iM8, iM8A, iM8B, iM9, iM9A, iM50, iM80; Goldway, Contec and Danatech, Penlon Intermed PM 8000, Vital Lite, Vital Line, Vital Wave.



sensor connection

< Cable length: 2,0m >




device connection

Mindray

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.MI2.IC2	

Compatible with: PM 5000, PM 6000, PM 8000.



sensor connection

< Cable length: 2,0m >



device connection

Mindray

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.MI3.IC2	





Compatible with: BeneView T5, T8 and iPM-8900.


Mindray

One piece ECG cable (3 or 5 leads)


Compatible with:

PM 7000, PM 8000, PM 9000 with 6 pin Connector, MEC 1000, MEC 1200.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



< Cable length: 3,5m >








Mindray

ECG trunk cable with leads (3 or 5 leads)


Compatible with:

PM 7000, PM 8000, PM 9000 with 6 pin Connector, MEC 1000, MEC 1200.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




Mindray

ECG trunk cable (3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	



< Cable length: 3,0m >







Compatible with: PM 7000, PM 8000, PM 9000 with 6 pin Connector, MEC 1000, MEC 1200.


Mindray

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
PM 7000, PM 8000, PM 9000
with 6 pin Connector, MEC
1000, MEC 1200.


Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >









Mindray

Temperature
probe.



< Cable length: 3,0m >



Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.1400PT	
Skin TITANIUM® 5,5mm		BIOTEMP.1400PTP	
Timpanic		BIOTEMP.1400T	
Adult esophageal / rectal		BIOTEMP.1400RA	
Pediatric esophageal / rectal		BIOTEMP.1400RP	


Compatible with:
T5 / T8 series, iPM-9800.






Mindray

Temperature
probe.



< Cable length: 3,0m >




Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.400PT	
Skin TITANIUM® 5,5mm		BIOTEMP.400PTP	
Timpanic		BIOTEMP.400T	
Adult esophageal / rectal		BIOTEMP.400RA	
Pediatric esophageal / rectal		BIOTEMP.400RP	

Compatible with:
All Mindray devices with
plug 1/4 connector.



< Cable length: 2,7m >




Moriya

SpO2 sensor.

Compatible with:

M1000, M3000, M9100.

Options		Code	
Adult Clip		BIOOX.MY.AC3	
Adult Silicone		BIOOX.MY.AS3	
Pediatric Clip		BIOOX.MY.PC3	
Pediatric Silicone		BIOOX.MY.PS3	
Neonate / Multisite		BIOOX.MY.Y3	
Wrap silicone		BIOOX.MY.WS3	



< Cable length: 0,9m >



Moriya


SpO2 short sensor.

this product requires an extension cable

Compatible with:

1001, 1002, 1003.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Moriya

SpO2 intermediary cable.

Compatible with: M1000, M3000, M9100.

Description	Code	
SpO2 intermediary cable	BIOOX.MO.IC2	


Moriya

One piece ECG cable
(3 or 5 leads)


Compatible with:

M2000, M3000, M9100.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	







< Cable length: 3,5m >




Moriya


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
M2000, M3000, M9100.

Options		Code
3 leads snap		BIOECG.BC.S3D
3 leads clamp		BIOECG.BC.C3D
5 leads snap		BIOECG.BC.S5D
5 leads clamp		BIOECG.BC.C5D






< Cable length: 4,0m >




Moriya

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code
3 leads		BIOECG.BC.3T
5 leads		BIOECG.BC.5T



< Cable length: 3,0m >







Compatible with: M2000, M3000, M9100.

Moriya

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
M2000, M3000, M9100.


Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5

Cable length: 1,0m >




Mortara

One piece EKG
10 lead cable.

Options		Code
Plug / Banana		BIOECG.MO.P10F




< Cable length: 3,0m >




Compatible with: ELI 10, ELI 50, ELI 100, ELI 150C, ELI 200, ELI 230, ELI 250C, ELI 350, Landscape, Portrait.


Mortara

One piece EKG
10 lead cable.


Options		Code
Plug / Banana		BIOECG.SC.P10F



< Cable length: 3,0m >



Compatible with: ELI 150, ELI 210, ELI 250.









< Cable length: 2,7m >




Nelcor


SpO2 sensor.

Compatible with:
N20, N180, N6000, NPB40, NPB75,
NPB3900, NPB190, NPB195, NPB4000.

Options		Code
Adult Clip		BIOOX.NC.AC3
Adult Silicone		BIOOX.NC.AS3
Pediatric Clip		BIOOX.NC.PC3
Pediatric Silicone		BIOOX.NC.PS3
Neonate / Multisite		BIOOX.NC.Y3
Wrap silicone		BIOOX.NC.WS3









< Cable length: 2,7m >




Nelcor


SpO2 sensor.

Compatible with:
N-100, N-200, N-1000 and N-2000
Preamp cable; MDE Escort E100,
Drager VitAlert 3200.

Options		Code
Adult Clip		BIOOX.NCL.AC3
Adult Silicone		BIOOX.NCL.AS3
Pediatric Clip		BIOOX.NCL.PC3
Pediatric Silicone		BIOOX.NCL.PS3
Neonate / Multisite		BIOOX.NCL.Y3
Wrap silicone		BIOOX.NCL.WS3



< Cable length: 2,7m >










Nellcor

SpO2 sensor.


Compatible with:

N-595 and N-600; CAS Med 511 Apnea Monitor; Cardell Veterinary 9403, 9500HD; Colin BP-S510 and T105NXM; Digicare LW600; Dinamap Carescape V100, ProCare, Pro Series, ProSeries V2, Pro 300, Pro 1000, Pro 1000 V3, E-NSAT Modules for i4 Modular Monitors with Oximax option; Fukuda Denshi DS-7100; GE Procure; Welch Allyn Spot Vital Signs 42NTB-E1, 42NT0-E1 LXi.

Options		Code	
Adult Clip		BIOOX.NCM.AC3	
Adult Silicone		BIOOX.NCM.AS3	
Pediatric Clip		BIOOX.NCM.PC3	
Pediatric Silicone		BIOOX.NCM.PS3	
Neonate / Multisite		BIOOX.NCM.Y3	
Wrap silicone		BIOOX.NCM.WS3	



< Cable length: 2,7m >









Nellcor

SpO2 sensor.

Compatible with:

Colin: BP88/306; Critikon Dinamap ProCare, Pro Series, Pro 300, Pro 1000, M-NSAT Modules; Datascope: Accutorr Plus, Duo, Passport 2, Passport XG, Spectrum, Trio; Nellcor: NPB-290, NPB-295, N-390, N-395, N-3000, N-6000, NPB-290, NPB-295, N-390, N-395, N-3000, N-6000; Nihon Kohden: BSM-2300, BSM-4100, BSM-5100, BSM- 4104A, BSM-4114A; Ohmeda enGuard CM5; Physio Control Lifepak 12.

Options		Code	
Adult Clip		BIOOX.NCN.AC3	
Adult Silicone		BIOOX.NCN.AS3	
Pediatric Clip		BIOOX.NCN.PC3	
Pediatric Silicone		BIOOX.NCN.PS3	
Neonate / Multisite		BIOOX.NCN.Y3	
Wrap silicone		BIOOX.NCN.WS3	



< Cable length: 2,7m >










Nellcor

SpO2 sensor.


Compatible with:

N 45, N 65, N 85, NPB 40, NPB 75, NBP-195; Colin; GE Dinamap Compact Dinamap MPS; All Welch Allyn . (with Oximax technology®)

Options		Code	
Adult Clip		BIOOX.NCO.AC3	
Adult Silicone		BIOOX.NCO.AS3	
Pediatric Clip		BIOOX.NCO.PC3	
Pediatric Silicone		BIOOX.NCO.PS3	
Neonate / Multisite		BIOOX.NCO.Y3	
Wrap silicone		BIOOX.NCO.WS3	









< Cable length: 0,9m >



Nellcor

SpO2 short sensor.
this product requires an extension cable

Compatible with:
All Nellcor® technology (standard) devices.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	









< Cable length: 0,9m >




Nellcor

SpO2 short sensor.
this product requires an extension cable

Compatible with:
All devices with Oximax technology®.

Options		Code	
Adult Clip		BIOOX.NCO.AC1	
Adult Silicone		BIOOX.NCO.AS1	
Pediatric Clip		BIOOX.NCO.PC1	
Pediatric Silicone		BIOOX.NCO.PS1	
Neonate / Multisite		BIOOX.NCO.Y1	
Wrap silicone		BIOOX.NCO.WS1	



sensor connection

< Cable length: 2,0m >




device connection

Nellcor

SpO2 intermediary cable.

Compatible with: N20, N180, N6000, NPB40, NPB75, NPB3900, NPB190, NPB195, NPB4000.

Description	Code	
SpO2 intermediary cable	BIOOX.NC.IC2	



sensor connection

< Cable length: 2,0m >




device connection

Nellcor

SpO2 intermediary cable.


Compatible with: N-100, N-200, N-1000 and N-2000 Preamp cable; MDE Escort E100, Drager VitAlert 3200.

Description	Code	
SpO2 intermediary cable	BIOOX.NCL.IC2	



sensor connection

< Cable length: 2,0m >




device connection

Nellcor

SpO2 intermediary cable.


Description	Code	
SpO2 intermediary cable	BIOOX.NCM.IC2	

Compatible with: N-595 and N-600; CAS Med 511 Apnea Monitor; Cardell Veterinary 9403, 9500HD; Colin BP-S510 and T105NXM; Digicare LW600; Dinamap Carescape V100, ProCare, Pro Series, ProSeries V2, Pro 300, Pro 1000, Pro 1000 V3, E-NSAT Modules for i4 Modular Monitors with Oximax option; Fukuda Denshi DS-7100; GE Procure; Welch Allyn Spot Vital Signs 42NTB-E1, 42NT0-E1 LXi.



sensor connection

< Cable length: 2,0m >




device connection

Nellcor

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.NCN.IC2	

Compatible with: Colin: BP88/306; Critikon Dinamap ProCare, Pro Series, Pro 300, Pro 1000, M-NSAT Modules; Datascope: Accutorr Plus, Duo, Passport 2, Passport XG, Spectrum, Trio; Nellcor: NPB-290, NPB-295, N-390, N-395, N-3000, N-6000, NPB-290, NPB-295, N-390, N-395, N-3000, N-6000; Nihon Kohden: BSM-2300, BSM-4100, BSM-5100, BSM- 4104A, BSM-4114A; Ohmeda enGuard CM5; Physio Control Lifepak 12.



sensor connection

< Cable length: 2,0m >



device connection

Nellcor

SpO2 intermediary cable.





Description	Code	
SpO2 intermediary cable	BIOOX.NCO.IC2	


Compatible with: N 45, N 65, N 85, NPB 40, NPB 75, NBP-195; Colin; GE Dinamap Compact Dinamap MPS; All Welch Allyn . (with Oximax technology®)

Nellcor


One piece ECG cable
(3 or 5 leads)

Compatible with:
N-200, N-250, C-Lock, N-1000, NPB-3900, NPB-4000.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	







< Cable length: 3,5m >




Nellcor


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
N-200, N-250, C-Lock, N-1000, NPB-3900, NPB-4000.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




Nellcor

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	



< Cable length: 3,0m >







Compatible with: N-200, N-250, C-Lock, N-1000, NPB-3900, NPB-4000.

Nellcor

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
N-200, N-250, C-Lock, N-1000, NPB-3900, NPB-4000.


Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >





< Cable length: 2,7m >










Nihon Kohden

SpO2 sensor.


Compatible with:

AL-800PA Module, BSM-1101 Life Scope EC, BSM-2101/2 Life Scope L, BSM-3101 Life Scope LC, BSM-7100/7200 Life Scope 8, BSM-8200/1 Life Scope 12, BSM-8301/2 Life Scope 9, SM-8800 Life Scope 14, OLV-1100A, NT3B.

Options		Code	
Adult Clip		BIOOX.NK.AC3	
Adult Silicone		BIOOX.NK.AS3	
Pediatric Clip		BIOOX.NK.PC3	
Pediatric Silicone		BIOOX.NK.PS3	
Neonate / Multisite		BIOOX.NK.Y3	
Wrap silicone		BIOOX.NK.WS3	



< Cable length: 2,7m >









Nihon Kohden

SpO2 sensor.

Compatible with:

BSM-4101, BSM-4102, BSM-4102A, 4103, 4111, 4113 Life Scope P, BSM-5105, 5135 Life Scope A, BSM-2351, 2353 Life Scope L, BSM-2301, 2303 Life Scope I, Life Scope J, BSM-9101 and PVM-2701.

Options		Code	
Adult Clip		BIOOX.NKT.AC3	
Adult Silicone		BIOOX.NKT.AS3	
Pediatric Clip		BIOOX.NKT.PC3	
Pediatric Silicone		BIOOX.NKT.PS3	
Neonate / Multisite		BIOOX.NKT.Y3	
Wrap silicone		BIOOX.NKT.WS3	



< Cable length: 0,9m >









Nihon Kohden


SpO2 short sensor.

this product requires an extension cable

Compatible with:

BSM-2300 Life Scope I, BSM-2301, MSM 2301K, 2303 Life Scope I, 2303K, BSM 2350 Life Scope L, BSM-2351, 2353 Life Scope L, BSM-4100 Life Scope P, BSM 4101, 4103 TL-101T, 4111, 4113 Life Scope P, BSM-5100 Life Scope A, BSM-5105, 5135 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM- 6501K, BSM-6701K, BSM-9510 Life Scope M, BSM-9800 Life Scope S, Oxipal P225F.

Options		Code	
Adult Clip		BIOOX.NKB.AC1	
Adult Silicone		BIOOX.NKB.AS1	
Pediatric Clip		BIOOX.NKB.PC1	
Pediatric Silicone		BIOOX.NKB.PS1	
Neonate / Multisite		BIOOX.NKB.Y1	
Wrap silicone		BIOOX.NKB.WS1	



sensor connection

< Cable length: 2,0m >




device connection

Nihon Kohden

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.NK.IC2	

Compatible with: AL-800PA Module, BSM-1101 Life Scope EC, BSM-2101/2 Life Scope L, BSM-3101 Life Scope LC, BSM-7100/7200 Life Scope 8, BSM-8200/1 Life Scope 12, BSM-8301/2 Life Scope 9, SM-8800 Life Scope 14, OLV-1100A, NT3B.



sensor connection

< Cable length: 2,0m >



device connection

Nihon Kohden

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.NKT.IC2	





Compatible with: BSM-4101, BSM-4102, BSM-4102A, 4103, 4111, 4113 Life Scope P, BSM-5105, 5135 Life Scope A, BSM-2351, 2353 Life Scope L, BSM-2301, 2303 Life Scope I, Life Scope J, BSM-9101 and PVM-2701.


Nihon Kohden

One piece ECG cable
(3 or 5 leads)


Compatible with:
Life Scope 6, Life Scope 8,

Life Scope 9, Life Scope 11, 6102, 6205, 6120, OEC 6102A (S/N < 20156), OEC-6105 (S/N < 20226), OEC-6120 A/J/K.

Options		Code
3 leads snap		BIOECG.NK.S3F
3 leads clamp		BIOECG.NK.C3F
5 leads snap		BIOECG.NK.S5F
5 leads clamp		BIOECG.NK.C5F



< Cable length: 3,5m >








Nihon Kohden

ECG trunk cable with leads
(3 or 5 leads)


Compatible with:
Life Scope 6, Life Scope 8,

Life Scope 9, Life Scope 11, 6102, 6205, 6120, OEC 6102A (S/N < 20156), OEC-6105 (S/N < 20226), OEC-6120 A/J/K.

Options		Code
3 leads snap		BIOECG.NK.S3D
3 leads clamp		BIOECG.NK.C3D
5 leads snap		BIOECG.NK.S5D
5 leads clamp		BIOECG.NK.C5D





< Cable length: 4,0m >




Nihon Kohden

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Compatible with: Life Scope 6, Life Scope 8, Life Scope 9, Life Scope 11, 6102, 6205, 6120, OEC 6102A (S/N < 20156), OEC-6105 (S/N < 20226), OEC-6120 A/J/K.

Options		Code
3 leads		BIOECG.NK.3T
5 leads		BIOECG.NK.5T







< Cable length: 3,0m >



Nihon Kohden

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
Life Scope 6, Life Scope 8, Life Scope 9, Life Scope 11, 6102, 6205, 6120, OEC 6102A (S/N < 20156), OEC-6105 (S/N < 20226), OEC-6120 A/J/K.

Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5

Cable length: 1,0m >




Nihon Kohden

One piece ECG cable
(3 or 5 leads)

Compatible with:

Merlin 78352C, Merlin 78354C, Merlin 78834C, M1001A, M1002A, Viridia, Omnicare CMA, ML65A, 66A, Heartstream and Codemaster XL, Tec-7621K, Tec-7631K, Tec-7721C, Tec-7731C, Tec-7700, Tec-5521, Tec-5531.




Options		Code
3 leads snap		BIOECG.NKV.S3F
3 leads clamp		BIOECG.NKV.C3F
5 leads snap		BIOECG.NKV.S5F
5 leads clamp		BIOECG.NKV.C5F

Nihon Kohden

ECG trunk cable with leads
(3 or 5 leads)

Compatible with:

Merlin 78352C, Merlin 78354C, Merlin 78834C, M1001A, M1002A, Viridia, Omnicare CMA, ML65A, 66A, Heartstream and Codemaster XL, Tec-7621K, Tec-7631K, Tec-7721C, Tec-7731C, Tec-7700, Tec-5521, Tec-5531.




Options		Code
3 leads snap		BIOECG.NKV.S3D
3 leads clamp		BIOECG.NKV.C3D
5 leads snap		BIOECG.NKV.S5D
5 leads clamp		BIOECG.NKV.C5D

Nihon Kohden

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Compatible with:

Merlin 78352C, Merlin 78354C, Merlin 78834C, M1001A, M1002A, Viridia, Omnicare CMA, ML65A, 66A, Heartstream and Codemaster XL, Tec-7621K, Tec-7631K, Tec-7721C, Tec-7731C, Tec-7700, Tec-5521, Tec-5531.



Options		Code
3 leads		BIOECG.NKV.3T
5 leads		BIOECG.NKV.5T


Nihon Kohden

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:

Merlin 78352C, Merlin 78354C,
Merlin 78834C, M1001A, M1002A, Viridia, Omnicare CMA, ML65A, 66A, Heartstream and Codemaster XL, Tec-7621K, Tec-7631K, Tec-7721C, Tec-7731C, Tec-7700, Tec-5521, Tec-5531.





Cable length: 1,0m >



Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5


Nihon Kohden

One piece ECG cable
(3 or 5 leads)


Options		Code
3 leads snap		BIOECG.NKT.S3F
3 leads clamp		BIOECG.NKT.C3F
5 leads snap		BIOECG.NKT.S5F
5 leads clamp		BIOECG.NKT.C5F

Compatible with:

BSM-2300 Life Scope I, BSM-2301A, BSM-2301K, BSM-2351A, BSM-4100, BSM-4102, BSM-4102A, Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701, BSM-9510 Life Scope M, BSM-9500 Life Scope S.







< Cable length: 3,5m >




Nihon Kohden

ECG trunk cable with leads
(3 or 5 leads)


Options		Code
3 leads snap		BIOECG.NKT.S3D
3 leads clamp		BIOECG.NKT.C3D
5 leads snap		BIOECG.NKT.S5D
5 leads clamp		BIOECG.NKT.C5D

Compatible with:

BSM-2300 Life Scope I, BSM-2301A, BSM-2301K, BSM-2351A, BSM-4100, BSM-4102, BSM-4102A, Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701, BSM-9510 Life Scope M, BSM-9500 Life Scope S.





< Cable length: 4,0m >




Nihon Kohden

ECG trunk cable
(3 or 5 leads).
this product requires a lead set


Options		Code
3 leads		BIOECG.NKT.3T
5 leads		BIOECG.NKT.5T

Compatible with:

BSM-2300 Life Scope I, BSM-2301A, BSM-2301K, BSM-2351A, BSM-4100, BSM-4102, BSM-4102A, Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701, BSM-9510 Life Scope M, BSM-9500 Life Scope S.



< Cable length: 3,0m >







Nihon Kohden


ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
BSM-2300 Life Scope I, BSM-2301A,

BSM-2301K, BSM-2351A, BSM-4100, BSM-4102, BSM-4102A, Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701, BSM-9510 Life Scope M, BSM-9500 Life Scope S.


Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Cable length: 1,0m >




Nihon Kohden

One piece EKG
10 lead cable.




< Cable length: 3,0m >




Compatible with: Cardiofax 651, 6581, 8110, 9010, 9110, 9620, Cardiofax Q, 1350, ECG-1250A, ECG-9620, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9130, ECG-9132, ECG-9110K, ECG-9130K, Cardioline (new line), Del Mar Reynolds Voyager 12; Newtech 1206; Carewell 1103, Carewell AR 1200ADV.

Nihon Kohden

One piece EKG
10 lead cable.



< Cable length: 3,0m >




Compatible with: Cardiofax/Q-Cardiofax, 6151, 6353, 6851, 8110, 8270, 8330, 8330A, 8350, 8370, 8420, 8423, 8450, 8453, 8820, 8823, 8830, 9020, 9022, 9130, 9132, Escape 12, Edan Smart ECG, Danatech, Dr. Lee 120B.

Nihon Kohden



< Cable length: 3,0m >



Temperature probe.


Compatible with:

All Nihon Kohden devices with plug 1/4 connector.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.400PT	
Skin TITANIUM® 5,5mm		BIOTEMP.400PTP	
Timpanic		BIOTEMP.400T	
Adult esophageal / rectal		BIOTEMP.400RA	
Pediatric esophageal / rectal		BIOTEMP.400RP	



< Cable length: 2,7m >




Nonin

SpO2 sensor.

Compatible with:

8500, 8550, 8600, 8700, 8800, 2500, 2500A, 9840, 9600, 9700, 2120, 7500, Respirationics 902M Plus, Lifesense and PureSAT.

Options		Code	
Adult Clip		BIOOX.NN.AC3	
Adult Silicone		BIOOX.NN.AS3	
Pediatric Clip		BIOOX.NN.PC3	
Pediatric Silicone		BIOOX.NN.PS3	
Neonate / Multisite		BIOOX.NN.Y3	
Wrap silicone		BIOOX.NN.WS3	



< Cable length: 0,9m >




Nonin

SpO2 short sensor.
this product requires an extension cable

Compatible with:

8500, 8550, 8600, 8700, 8800, 2500, 2500A, 9840, 9600, 9700, 2120, 7500, Respirationics 902M Plus, Lifesense and PureSAT.

Options		Code	
Adult Clip		BIOOX.NN.AC1	
Adult Silicone		BIOOX.NN.AS1	
Pediatric Clip		BIOOX.NN.PC1	
Pediatric Silicone		BIOOX.NN.PS1	
Neonate / Multisite		BIOOX.NN.Y1	
Wrap silicone		BIOOX.NN.WS1	



sensor connection

< Cable length: 2,0m >




device connection

Nonin

SpO2 intermediary cable.


Description	Code	
SpO2 intermediary cable	BIOOX.NN.IC2	

Compatible with: 8500, 8550, 8600, 8700, 8800, 2500, 2500A, 9840, 9600, 9700, 2120, 7500, Respironics 902M Plus, Lifesense and PureSAT.



sensor connection

< Cable length: 2,0m >




device connection

Nonin

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.NN2.IC2	83.80

Compatible with: 8600, 8500, 8700, 8800, 96, 009, 700.



< Cable length: 2,7m >



Novamatrix

SpO2 sensor.

Compatible with:
505, 510, 511, 515, 515A, 515B, 515C, 520, 520A, 7000, 7100, Oxypleth; Spacelabs: 90351-D/G, 90369, 9046, 90466, 90467, 90489, 90651-07/08, 90419-01, 90651A-07.

Options		Code	
Adult Clip		BIOOX.SP.AC3	
Adult Silicone		BIOOX.SP.AS3	
Pediatric Clip		BIOOX.SP.PC3	
Pediatric Silicone		BIOOX.SP.PS3	
Neonate / Multisite		BIOOX.SP.Y3	
Wrap silicone		BIOOX.SP.WS3	



< Cable length: 0,9m >










Novamatrix

SpO2 short sensor.
this product requires an extension cable

Compatible with:

Respironics 920, 920M, 920M Plus, Tidal Wave 710Sp, 715Sp. Please note that some Respironics 920 models require a 7-pin connector, for these models please request item S403-08T0.

Options		Code	
Adult Clip		BIOOX.NOD.AC1	
Adult Silicone		BIOOX.NOD.AS1	
Pediatric Clip		BIOOX.NOD.PC1	
Pediatric Silicone		BIOOX.NOD.PS1	
Neonate / Multisite		BIOOX.NOD.Y1	
Wrap silicone		BIOOX.NOD.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Novamatrix


SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.NO.IC2	

Compatible with: 505, 510, 511, 515, 515A, 515B, 515C, 520, 520A, 7000, 7100, Oxypleth; Spacelabs: 90351-D/G, 90369, 9046, 90466, 90467, 90489, 90651-07/08, 90419-01, 90651A-07.


Ohio

Temperature probe.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.100PT	
Skin TITANIUM® 5,5mm		BIOTEMP.100PTP	




< Cable length: 1,5m >



Compatible with: 190.



< Cable length: 2,7m >



Pace Tech

SpO2 sensor.

Compatible with:

Generra: 710, 200, 300, 400, 810, 530;
Pace Tech: 520.

Options		Code	
Adult Clip		BIOOX.PT.AC3	
Adult Silicone		BIOOX.PT.AS3	
Pediatric Clip		BIOOX.PT.PC3	
Pediatric Silicone		BIOOX.PT.PS3	
Neonate / Multisite		BIOOX.PT.Y3	
Wrap silicone		BIOOX.PT.WS3	



< Cable length: 2,7m >



Pace Tech

SpO2 sensor.

Compatible with:

Vitalmax 4000, Vitalmax 4100 (1st version).

Options		Code	
Adult Clip		BIOOX.PTR.AC3	
Adult Silicone		BIOOX.PTR.AS3	
Pediatric Clip		BIOOX.PTR.PC3	
Pediatric Silicone		BIOOX.PTR.PS3	
Neonate / Multisite		BIOOX.PTR.Y3	
Wrap silicone		BIOOX.PTR.WS3	



< Cable length: 2,7m >




Pace Tech

SpO2 sensor.


Compatible with:

Mini Pack 300.

Options		Code	
Adult Clip		BIOOX.BC.AC3	
Adult Silicone		BIOOX.BC.AS3	
Pediatric Clip		BIOOX.BC.PC3	
Pediatric Silicone		BIOOX.BC.PS3	
Neonate / Multisite		BIOOX.BC.Y3	
Wrap silicone		BIOOX.BC.WS3	



< Cable length: 0,9m >



Pace Tech


SpO2 short sensor.

this product requires an extension cable

Compatible with:

Generra: 710, 200, 300, 400, 810, 530; Pace Tech: 520, Vitalmax 4000, Vitalmax 4100 (1st version).

Options		Code	
Adult Clip		BIOOX.BC.AC1	
Adult Silicone		BIOOX.BC.AS1	
Pediatric Clip		BIOOX.BC.PC1	
Pediatric Silicone		BIOOX.BC.PS1	
Neonate / Multisite		BIOOX.BC.Y1	
Wrap silicone		BIOOX.BC.WS1	



sensor connection

< Cable length: 2,0m >




device connection

Pace Tech

SpO2 intermediary cable.

Compatible with: Generra: 710, 200, 300, 400, 810, 530; Pace Tech: 520.

Description	Code	
SpO2 intermediary cable	BIOOX.PT.IC2	



sensor connection

< Cable length: 2,0m >



device connection

Pace Tech


SpO2 intermediary cable.

Compatible with: Vitalmax 4000, Vitalmax 4100 (1st version).

Description	Code	
SpO2 intermediary cable	BIOOX.PTR.IC2	



< Cable length: 2,7m >




Palco

SpO2 sensor.


Compatible with:

All Palco devices.

Options		Code	
Adult Clip		BIOOX.PA.AC3	
Adult Silicone		BIOOX.PA.AS3	
Pediatric Clip		BIOOX.PA.PC3	
Pediatric Silicone		BIOOX.PA.PS3	
Neonate / Multisite		BIOOX.PA.Y3	
Wrap silicone		BIOOX.PA.WS3	



< Cable length: 2,7m >










Philips - HP

SpO2 sensor.


Compatible with:

78352A, 78352C, 78354A, 78354C, 78834C, M1020A, M1025A, M1025B, M1092A, M1094B, M1204A, M1205A, M1350B, M1722A, M1722B, M1190A, M1723A, M1723B, M1732A, M1732B, M2475B, M3046A M4, V24E, Codemaster, HeartStart XL, Agilent 50XM Fetal Monitor.

Options		Code	
Adult Clip		BIOOX.PH1.AC3	
Adult Silicone		BIOOX.PH1.AS3	
Pediatric Clip		BIOOX.PH1.PC3	
Pediatric Silicone		BIOOX.PH1.PS3	
Neonate / Multisite		BIOOX.PH1.Y3	
Wrap silicone		BIOOX.PH1.WS3	



< Cable length: 2,7m >









Philips - HP

SpO2 sensor.

Compatible with:

C1, C3, M3, M4, MP20, MP30, MP40, MP50, MP60, MP70, MP80, MP90, VM4, VM6, VM8, HP Viridia 24C, Philips HeartStart defibrillator and IntelliVue; Goldway G40; Solaris NT2A, NT2C.

Options		Code	
Adult Clip		BIOOX.PH2.AC3	
Adult Silicone		BIOOX.PH2.AS3	
Pediatric Clip		BIOOX.PH2.PC3	
Pediatric Silicone		BIOOX.PH2.PS3	
Neonate / Multisite		BIOOX.PH2.Y3	
Wrap silicone		BIOOX.PH2.WS3	



< Cable length: 2,7m >










Philips - HP

SpO2 sensor.


Compatible with:

MP20, MP30, MP40, MP50, MP60, MP70, MP90, VN4, VN6, VN8 with M. Tech.

Options		Code	
Adult Clip		BIOOX.PH3.AC3	
Adult Silicone		BIOOX.PH3.AS3	
Pediatric Clip		BIOOX.PH3.PC3	
Pediatric Silicone		BIOOX.PH3.PS3	
Neonate / Multisite		BIOOX.PH3.Y3	
Wrap silicone		BIOOX.PH3.WS3	



< Cable length: 0,9m >



Philips - HP

SpO2 short sensor.

this product requires an extension cable

Compatible with:

Viridia 24C, C1, C3, M3, M4, MP20, MP30, MP40, MP50, MP60, MP70, MP80, MP90, VM4, VM6, VM8, HeartStart MRx and SureSigns VS3.

Options		Code	
Adult Clip		BIOOX.PH2.AC1	
Adult Silicone		BIOOX.PH2.AS1	
Pediatric Clip		BIOOX.PH2.PC1	
Pediatric Silicone		BIOOX.PH2.PS1	
Neonate / Multisite		BIOOX.PH2.Y1	
Wrap silicone		BIOOX.PH2.WS1	



< Cable length: 0,9m >



Philips - HP


SpO2 short sensor.

this product requires an extension cable

Compatible with:

78352A, 78352C, 78354A, 78354C, 78834C, M1020A, M1025A, M1025B, M1092A, M1094B, M1204A, M1205A, M1350B, M1722A, M1722B, M1190A, M1723A, M1723B, M1732A, M1732B, M2475B, M3046A M4, V24E, Codemaster, HeartStart XL, Agilent 50XM Fetal Monitor, C1, C3, M3, M4, MP20, MP30, MP40, MP50, MP60, MP70, MP80, MP90, VM4, VM6, VM8, HP Viridia 24C, Philips HeartStart defibrillator and IntelliVue; Goldway G40; Solaris NT2A, NT2C.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >




device connection

Philips - HP

SpO2 intermediary cable.


Description	Code	
SpO2 intermediary cable	BIOOX.PH1.IC2	

Compatible with: 78352A, 78352C, 78354A, 78354C, 78834C, M1020A, M1025A, M1025B, M1092A, M1094B, M1204A, M1205A, M1350B, M1722A, M1722B, M1190A, M1723A, M1723B, M1732A, M1732B, M2475B, M3046A M4, V24E, Codemaster, HeartStart XL, Agilent 50XM Fetal Monitor.



sensor connection

< Cable length: 2,0m >



device connection

Philips - HP

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.PH2.IC2	

Compatible with: Viridia 24C, C1, C3, M3, M4, MP20, MP30, MP40, MP50, MP60, MP70, MP80, MP90, VM4, VM6, VM8, HeartStart MRx and SureSigns VS3.


Philips - HP

One piece ECG cable (3 or 5 leads)


Options		Code	
3 leads snap		BIOECG.HP1.S3F	
3 leads clamp		BIOECG.HP1.C3F	
5 leads snap		BIOECG.HP1.S5F	
5 leads clamp		BIOECG.HP1.C5F	

Compatible with:

Codemaster XL, Codemaster 100, 43100A, 43110A, 43110M, 43120A, MC Defibrillator, 43200A Portable, M3500B, 77000AC, 77020A, 77010AC, Sonos/AC/AR, 77025A, 77030A, 77035A, 77040A, 78208A, 78210A, 78330A, 7833 1A, 78332A, 78333A, 78341A, 78342A, 78345A, 78346A, 78351T, 78352A/C, 78353A/B, 78354A, 78532A/B, 78534A/B/C, 78619A, 78620A, 78660A/B, 78661A/B, 78670A, 78671A, M1722A/B, M1723A/B, M2400A, M2406A, Sonos 2000/2500, tems, M2425A, M2475B.




< Cable length: 3,5m >




Philips - HP

ECG trunk cable with leads (3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.HP1.S3D	
3 leads clamp		BIOECG.HP1.C3D	
5 leads snap		BIOECG.HP1.S5D	
5 leads clamp		BIOECG.HP1.C5D	



< Cable length: 4,0m >






Compatible with:

Codemaster XL, Codemaster 100, 43100A, 43110A, 43110M, 43120A, MC Defibrillator, 43200A Portable, M3500B, 77000AC, 77020A, 77010AC, Sonos/AC/AR, 77025A, 77030A, 77035A, 77040A, 78208A, 78210A, 78330A, 7833 1A, 78332A, 78333A, 78341A, 78342A, 78345A, 78346A, 78351T, 78352A/C, 78353A/B, 78354A, 78532A/B, 78534A/B/C, 78619A, 78620A, 78660A/B, 78661A/B, 78670A, 78671A, M1722A/B, M1723A/B, M2400A, M2406A, Sonos 2000/2500, tems, M2425A, M2475B.


Philips - HP

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.HP1.3T	
5 leads		BIOECG.HP1.5T	



< Cable length: 3,0m >







Compatible with:

Codemaster XL, Codemaster 100, 43100A, 43110A, 43110M, 43120A, MC Defibrillator, 43200A Portable, M3500B, 77000AC, 77020A, 77010AC, Sonos/AC/AR, 77025A, 77030A, 77035A, 77040A, 78208A, 78210A, 78330A, 7833 1A, 78332A, 78333A, 78341A, 78342A, 78345A, 78346A, 78351T, 78352A/C, 78353A/B, 78354A, 78532A/B, 78534A/B/C, 78619A, 78620A, 78660A/B, 78661A/B, 78670A, 78671A, M1722A/B, M1723A/B, M2400A, M2406A, Sonos 2000/2500, tems, M2425A, M2475B.

Philips - HP


ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Compatible with:





Codemaster XL, Codemaster 100, 43100A, 43110A, 43110M, 43120A, MC Defibrillator, 43200A Portable, M3500B, 77000AC, 77020A, 77010AC, Sonos/AC/AR, 77025A, 77030A, 77035A, 77040A, 78208A, 78210A, 78330A, 7833 1A, 78332A, 78333A, 78341A, 78342A, 78345A, 78346A, 78351T, 78352A/C, 78353A/B, 78354A, 78532A/B, 78534A/B/C, 78619A, 78620A, 78660A/B, 78661A/B, 78670A, 78671A, M1722A/B, M1723A/B, M2400A, M2406A, Sonos 2000/2500, tems, M2425A, M2475B.

Cable length: 1,0m >




Philips - HP

One piece ECG cable
(3 or 5 leads)


Options		Code	
3 leads snap		BIOECG.HP2.S3F	
3 leads clamp		BIOECG.HP2.C3F	
5 leads snap		BIOECG.HP2.S5F	
5 leads clamp		BIOECG.HP2.C5F	

Compatible with:

Merlin, Viridia, Omnicare CMA, Heartstream, Codemaster XL, Codemaster 100, Philips Heartstart desfibrilator, Philips Agilent, Envisor, IntelliVue, MP-20, MP-30, MP-50, IE33, 43100A, 43110A, 43120A, 78352C, 78354C, 78834C, 862474-C3, M1001A/B, M1002A/B, M3002A, MMS X2, M8102A MP2, M8105A MP5, M1165A, M1166A, M1167A, M1175A, M1177A, M1204A, M1205A, M1264A, M1275A Transport, M1280A, M1722A/B, M1723A/B, M2400A, M2406A, M2601A, M3000A, M3001A, M3046A, M2, M3, M4, M4735A, V24E, Agilent 50XM Fetal Monitor, Newtech NT3A, 862472, 862473, 862474, 862475, 862478, 862479, 863064, 863065, 863066, 863068, VM6, VM8, M4735A, HeartStart XL, M3535A, M3536A, Solaris NT3A, NT3B, NT3C, NT3E, NT3F.







< Cable length: 3,5m >




Philips - HP

ECG trunk cable with leads
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.HP2.S3D	
3 leads clamp		BIOECG.HP2.C3D	
5 leads snap		BIOECG.HP2.S5D	
5 leads clamp		BIOECG.HP2.C5D	

Compatible with:



Merlin, Viridia, Omnicare CMA, Heartstream, Codemaster XL, Codemaster 100, Philips Heartstart defibrillator, Philips Agilent, Envisor, IntelliVue, MP-20, MP-30, MP-50, IE33, 43100A, 43110A, 43120A, 78352C, 78354C, 78834C, 862474-C3, M1001A/B, M1002A/B, M3002A, MMS X2, M8102A MP2, M8105A MP5, M1165A, M1166A, M1167A, M1175A, M1177A, M1204A, M1205A, M1264A, M1275A Transport, M1280A, M1722A/B, M1723A/B, M2400A, M2406A, M2601A, M3000A, M3001A, M3046A, M2, M3, M4, M4735A, V24E, Agilent 50XM Fetal Monitor, Newtech NT3A, 862472, 862473, 862474, 862475, 862478, 862479, 863064, 863065, 863066, 863068, VM6, VM8, M4735A, HeartStart XL, M3535A, M3536A, Solaris NT3A, NT3B, NT3C, NT3E, NT3F.



Philips - HP

ECG trunk cable
(3 or 5 leads).

this product requires a lead set

Options		Code	
3 leads		BIOECG.HP2.3T	
5 leads		BIOECG.HP2.5T	







Compatible with:

Merlin, Viridia, Omnicare CMA, Heartstream, Codemaster XL, Codemaster 100, Philips Heartstart defibrillator, Philips Agilent, Envisor, IntelliVue, MP-20, MP-30, MP-50, IE33, 43100A, 43110A, 43120A, 78352C, 78354C, 78834C, 862474-C3, M1001A/B, M1002A/B, M3002A, MMS X2, M8102A MP2, M8105A MP5, M1165A, M1166A, M1167A, M1175A, M1177A, M1204A, M1205A, M1264A, M1275A Transport, M1280A, M1722A/B, M1723A/B, M2400A, M2406A, M2601A, M3000A, M3001A, M3046A, M2, M3, M4, M4735A, V24E, Agilent 50XM Fetal Monitor, Newtech NT3A, 862472, 862473, 862474, 862475, 862478, 862479, 863064, 863065, 863066, 863068, VM6, VM8, M4735A, HeartStart XL, M3535A, M3536A, Solaris NT3A, NT3B, NT3C, NT3E, NT3F.

Philips - HP


ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Compatible with:


Merlin, Viridia, Omnicare CMA, Heartstream, Codemaster XL, Codemaster 100, Philips Heartstart desfibrilator, Philips Agilent, Envisor, IntelliVue, MP-20, MP-30, MP-50, IE33, 43100A, 43110A, 43120A, 78352C, 78354C, 78834C, 862474-C3, M1001A/B, M1002A/B, M3002A, MMS X2, M8102A MP2, M8105A MP5, M1165A, M1166A, M1167A, M1175A, M1177A, M1204A, M1205A, M1264A, M1275A Transport, M1280A, M1722A/B, M1723A/B, M2400A, M2406A, M2601A, M3000A, M3001A, M3046A, M2, M3, M4, M4735A, V24E, Agilent 50XM Fetal Monitor, Newtech NT3A, 862472, 862473, 862474, 862475, 862478, 862479, 863064, 863065, 863066, 863068, VM6, VM8, M4735A, HeartStart XL, M3535A, M3536A, Solaris NT3A, NT3B, NT3C, NT3E, NT3F.


Cable length: 1,0m >




Philips - HP

One piece EKG
10 lead cable.

Options		Code	
Plug / Banana		BIOECG.HP.P10F	139.80




< Cable length: 3,0m >



Compatible with: M1770A, M1771A, M1772A, PageWriter 100, 200, 200i and Agilent 300i, 300PI; CardioTech CRT-6U, CRT-6D, CRT-6T; Bionet Cardiocare 2000; Delmar Reynolds CardioCollect 12; BLT-08 LT ECG; Spacelabs CardioDirect 12; Woodway Cortex CPX System; Welch Allyn CP 100 and CP 200.








< Cable length: 3,0m >



Philips - HP

Temperature probe.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.1200PT	
Skin TITANIUM® 5,5mm		BIOTEMP.1200PTP	
Timpanic		BIOTEMP.1200T	
Adult esophageal / rectal		BIOTEMP.1200RA	
Pediatric esophageal / rectal		BIOTEMP.1200RP	

Compatible with:





Módulo de probe M1029.


Physio Control / Medtronic

One piece ECG cable
(3 or 5 leads)


Compatible with:

Lifepak 5, Lifepak 6, Lifepak 6S, Lifepak 7, Lifepak 8, Lifepak 9, Lifepak 9P, Lifepak 10, VSM-1, VSM-2, VSM-2ESF, VSM-3, 1440.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



< Cable length: 3,5m >







Physio Control / Medtronic

ECG trunk cable with leads
(3 or 5 leads)


Compatible with:

Lifepak 5, Lifepak 6, Lifepak 6S, Lifepak 7, Lifepak 8, Lifepak 9, Lifepak 9P, Lifepak 10, VSM-1, VSM-2, VSM-2ESF, VSM-3, 1440.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




Physio Control / Medtronic

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	



< Cable length: 3,0m >



Compatible with:





Lifepak 5, Lifepak 6, Lifepak 6S, Lifepak 7, Lifepak 8, Lifepak 9, Lifepak 9P, Lifepak 10, VSM-1, VSM-2, VSM-2ESF, VSM-3, 1440.

Physio Control / Medtronic

ECG lead set


(3 or 5 leads)
this product requires a trunk cable

Compatible with:

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	





Lifepak 5, Lifepak 6, Lifepak 6S, Lifepak 7, Lifepak 8, Lifepak 9, Lifepak 9P, Lifepak 10, VSM-1, VSM-2, VSM-2ESF, VSM-3, 1440.


Cable length: 1,0m >




Physio Control / Medtronic

One piece ECG cable
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.PC.S3F	
3 leads clamp		BIOECG.PC.C3F	
5 leads snap		BIOECG.PC.S5F	
5 leads clamp		BIOECG.PC.C5F	







< Cable length: 3,5m >




Compatible with: Lifepak 11, Lifepak 12, Lifepak 15, Lifepak 20.


Physio Control / Medtronic

ECG trunk cable with leads
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.PC.S3D	
3 leads clamp		BIOECG.PC.C3D	
5 leads snap		BIOECG.PC.S5D	
5 leads clamp		BIOECG.PC.C5D	



< Cable length: 4,0m >






Compatible with: Lifepak 11, Lifepak 12, Lifepak 15, Lifepak 20.

Physio Control / Medtronic


ECG trunk cable

(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.PC.3T	
5 leads		BIOECG.PC.5T	



< Cable length: 3,0m >







Compatible with: Lifepak 11, Lifepak 12, Lifepak 15, Lifepak 20.

Physio Control / Medtronic

ECG lead set


(3 or 5 leads)
this product requires a trunk cable

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Compatible with:


Lifepak 11, Lifepak 12, Lifepak 15,
Lifepak 20.


Cable length: 1,0m >




Progetti

One piece EKG
10 lead cable.


Options		Code	
Plug / Banana		BIOECG.SC2.P10F	




< Cable length: 3,0m >



Compatible with: EPG Smart, EPG Pocket, EPG Junior, EPG View, EPG 6 View, EPG Plus, EPG 612, EPG 13.



< Cable length: 2,7m >









Protec

SpO2 sensor.


Compatible with:

Provision 4000, Provision 7000.

Options		Code	
Adult Clip		BIOOX.PR.AC3	
Adult Silicone		BIOOX.PR.AS3	
Pediatric Clip		BIOOX.PR.PC3	
Pediatric Silicone		BIOOX.PR.PS3	
Neonate / Multisite		BIOOX.PR.Y3	
Wrap silicone		BIOOX.PR.WS3	









< Cable length: 0,9m >




Protec

SpO2 short sensor.
this product requires an extension cable

Compatible with:
Provision 4000, Provision 7000.

Options		Code	
Adult Clip		BIOOX.BC.AC1	
Adult Silicone		BIOOX.BC.AS1	
Pediatric Clip		BIOOX.BC.PC1	
Pediatric Silicone		BIOOX.BC.PS1	
Neonate / Multisite		BIOOX.BC.Y1	
Wrap silicone		BIOOX.BC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Protec


SpO2 intermediary cable.

Compatible with: Provision 4000, Provision 7000.

Description	Code	
SpO2 intermediary cable	BIOOX.PR.IC2	





< Cable length: 1,5m >



Rosch


Temperature probe.

Compatible with: All Rosch devices.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.850PT	59.80
Skin TITANIUM® 5,5mm		BIOTEMP.850PTP	59.80


Schiller

One piece EKG 10 lead cable.

Options		Code	
Plug / Banana		BIOECG.SC.P10F	




< Cable length: 3,0m >




Compatible with: AT3, AT6, CS6, AT5, AT10, AT60, CS-200; Bionet Cardiotouch.


Schiller

One piece EKG
10 lead cable.

Options		Code
Plug / Banana		BIOECG.SC2.P10F




< Cable length: 3,0m >



Compatible with: AT-1, AT-2, AT-2 Light, AT-2 Plus, AT-4, AT-10, AT-10 Plus, AT-60, AT-101, AT-102, AT-102-Plus, AT-104, AT-110, CS-100, CS-200, MS-2010, MS-2015; Bionet: CardioCare 2000; Esaote: P80 Series 40893, Biosound MyLab 30 CD, 30 CV; Cardiette: AR1200 adv; Welch Allyn CP10, CP20.









< Cable length: 2,7m >



Spacelabs


SpO2 sensor.

Compatible with:
90351-D/G, 90369, 9046, 90466, 90467,
90489, 90651-07/08, 90419-01, 90651A-07.

Options		Code
Adult Clip		BIOOX.SP.AC3
Adult Silicone		BIOOX.SP.AS3
Pediatric Clip		BIOOX.SP.PC3
Pediatric Silicone		BIOOX.SP.PS3
Neonate / Multisite		BIOOX.SP.Y3
Wrap silicone		BIOOX.SP.WS3









< Cable length: 2,7m >




Spacelabs

SpO2 sensor.

Compatible with:
90303B, 90351-0/6, 90465, 0466, 90489,
90651A-08, 90651-A, IM77, 90308, Scout.

Options		Code
Adult Clip		BIOOX.SPN.AC3
Adult Silicone		BIOOX.SPN.AS3
Pediatric Clip		BIOOX.SPN.PC3
Pediatric Silicone		BIOOX.SPN.PS3
Neonate / Multisite		BIOOX.SPN.Y3
Wrap silicone		BIOOX.SPN.WS3



< Cable length: 2,7m >




Spacelabs

SpO2 sensor.

Compatible with:

1050, 1600, 1700, 90367, 90369, 90469, 90496, 90651-08, Ultraview, 91220, 91369, mCare 300, sólo con tecnología Nellcor®.

Options		Code	
Adult Clip		BIOOX.SPR.AC3	
Adult Silicone		BIOOX.SPR.AS3	
Pediatric Clip		BIOOX.SPR.PC3	
Pediatric Silicone		BIOOX.SPR.PS3	
Neonate / Multisite		BIOOX.SPR.Y3	
Wrap silicone		BIOOX.SPR.WS3	



< Cable length: 0,9m >



Spacelabs


SpO2 short sensor.

this product requires an extension cable

Compatible with:

All Spacelabs devices with Nellcor® technology (standard).

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >




device connection

Spacelabs

SpO2 intermediary cable.

Description	Code	
SpO2 intermediary cable	BIOOX.SPN.IC2	

Compatible with: 90303B, 90351-0/6, 90465, 0466, 90489, 90651A-08, 90651-A, IM77, 90308, Scout.



sensor connection

< Cable length: 2,0m >



device connection

Spacelabs





SpO2 intermediary cable.


Description	Code	
SpO2 intermediary cable	BIOOX.SPR.IC2	

Compatible with: 1050, 1600, 1700, 90367, 90369, 90469, 90496, 90651-08, Ultraview, 91220, 91369, mCare 300, Nellcor® technology ONLY.


Spacelabs

One piece ECG cable
(3 or 5 leads)

Options		Code
3 leads snap		BIOECG.BC.S3F
3 leads clamp		BIOECG.BC.C3F
5 leads snap		BIOECG.BC.S5F
5 leads clamp		BIOECG.BC.C5F







< Cable length: 3,5m >




Compatible with: 500 / 600 / 900 Series (w/o Respiration), 90300 / 90400 / 90600 / 90700 Series, Tek-208 / 408 / 511 / 512 514 / 521, PC Bedsi, PCMS.


Spacelabs

ECG trunk cable with leads
(3 or 5 leads)

Options		Code
3 leads snap		BIOECG.BC.S3D
3 leads clamp		BIOECG.BC.C3D
5 leads snap		BIOECG.BC.S5D
5 leads clamp		BIOECG.BC.C5D





< Cable length: 4,0m >




Compatible with: 500 / 600 / 900 Series (w/o Respiration), 90300 / 90400 / 90600 / 90700 Series, Tek-208 / 408 / 511 / 512 514 / 521, PC Bedsi, PCMS.


Spacelabs

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code
3 leads		BIOECG.BC.3T
5 leads		BIOECG.BC.5T



< Cable length: 3,0m >







Compatible with: 500 / 600 / 900 Series (w/o Respiration), 90300 / 90400 / 90600 / 90700 Series, Tek-208 / 408 / 511 / 512 514 / 521, PC Bedsi, PCMS.

Spacelabs

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
500 / 600 / 900 Series (w/o Respiration), 90300 / 90400 / 90600 / 90700 Series, Tek-208 / 408 / 511 / 512 514 / 521, PC Bedsi, PCMS.





Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	


Cable length: 1,0m >




Spacelabs

One piece ECG cable
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.SP.S3F	
3 leads clamp		BIOECG.SP.C3F	
5 leads snap		BIOECG.SP.S5F	
5 leads clamp		BIOECG.SP.C5F	







< Cable length: 3,5m >




Compatible with: Ultraview 90309, 90367, 90492, 90496, 91220, 91496.


Spacelabs

ECG trunk cable with leads
(3 or 5 leads)

Options		Code	
3 leads snap		BIOECG.SP.S3D	
3 leads clamp		BIOECG.SP.C3D	
5 leads snap		BIOECG.SP.S5D	
5 leads clamp		BIOECG.SP.C5D	





< Cable length: 4,0m >




Compatible with: Ultraview 90309, 90367, 90492, 90496, 91220, 91496.

Spacelabs

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code
3 leads		BIOECG.SP.3T
5 leads		BIOECG.SP.5T







< Cable length: 3,0m >



Compatible with: Ultraview 90309, 90367, 90492, 90496, 91220, 91496.


Spacelabs

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5

Compatible with:
Ultraview 90309, 90367, 90492,
90496, 91220, 91496.

Cable length: 1,0m >









< Cable length: 2,7m >




Takaoka

SpO2 sensor.

Compatible with:
All Takaoka devices with
Redel 6 pin connector.

Options		Code
Adult Clip		BIOOX.TA.AC3
Adult Silicone		BIOOX.TA.AS3
Pediatric Clip		BIOOX.TA.PC3
Pediatric Silicone		BIOOX.TA.PS3
Neonate / Multisite		BIOOX.TA.Y3
Wrap silicone		BIOOX.TA.WS3









< Cable length: 0,9m >



Takaoka

SpO2 short sensor.
this product requires an extension cable

Compatible with:
All Takaoka devices.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

Takaoka

SpO2 intermediary cable.





Compatible with: All Takaoka devices with Redel 6 pin connector.


Description	Code	
SpO2 intermediary cable	BIOOX.TA.IC2	

Takaoka


One piece ECG cable
(3 or 5 leads)

Compatible with:
9100, 9200, Cardiotec, Oxifast.

Options		Code
3 leads snap		BIOECG.BC.S3F
3 leads clamp		BIOECG.BC.C3F
5 leads snap		BIOECG.BC.S5F
5 leads clamp		BIOECG.BC.C5F







< Cable length: 3,5m >




Takaoka


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
9100, 9200, Cardiotec, Oxifast.

Options		Code
3 leads snap		BIOECG.BC.S3D
3 leads clamp		BIOECG.BC.C3D
5 leads snap		BIOECG.BC.S5D
5 leads clamp		BIOECG.BC.C5D






< Cable length: 4,0m >




Takaoka

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code
3 leads		BIOECG.BC.3T
5 leads		BIOECG.BC.5T



< Cable length: 3,0m >







Compatible with: 9100, 9200, Cardiotec, Oxifast.

Takaoka


ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Compatible with:
9100, 9200, Cardiotec, Oxifast.

Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5

Cable length: 1,0m >





< Cable length: 2,7m >




TEB

SpO2 sensor.

Compatible with:
MPC 10, MPC 20.

Options		Code	
Adult Clip		BIOOX.TE.AC3	
Adult Silicone		BIOOX.TE.AS3	
Pediatric Clip		BIOOX.TE.PC3	
Pediatric Silicone		BIOOX.TE.PS3	
Neonate / Multisite		BIOOX.TE.Y3	
Wrap silicone		BIOOX.TE.WS3	



< Cable length: 0,9m >




TEB

SpO2 short sensor.

this product requires an extension cable

Compatible with:
MPC 10, MPC 20.

Options		Code	
Adult Clip		BIOOX.NC.AC1	
Adult Silicone		BIOOX.NC.AS1	
Pediatric Clip		BIOOX.NC.PC1	
Pediatric Silicone		BIOOX.NC.PS1	
Neonate / Multisite		BIOOX.NC.Y1	
Wrap silicone		BIOOX.NC.WS1	



sensor connection

< Cable length: 2,0m >



device connection

TEB

SpO2 intermediary cable.





Compatible with: MPC 10, MPC 20.


Description	Code	
SpO2 intermediary cable	BIOOX.TE.IC2	

TEB


One piece ECG cable
(3 or 5 leads)

Compatible with:
M10, M10B, Med1, M10A, M10C.

Options		Code
3 leads snap		BIOECG.TE.S3F
3 leads clamp		BIOECG.TE.C3F
5 leads snap		BIOECG.TE.S5F
5 leads clamp		BIOECG.TE.C5F







< Cable length: 3,5m >




TEB


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
M10, M10B, Med1, M10A, M10C.

Options		Code
3 leads snap		BIOECG.TE.S3D
3 leads clamp		BIOECG.TE.C3D
5 leads snap		BIOECG.TE.S5D
5 leads clamp		BIOECG.TE.C5D






< Cable length: 4,0m >




TEB

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code
3 leads		BIOECG.TE.3T
5 leads		BIOECG.TE.5T







< Cable length: 3,0m >



Compatible with: M10, M10B, Med1, M10A, M10C.

TEB

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5


Compatible with:
M10, M10B, Med1, M10A, M10C.

Cable length: 1,0m >





< Cable length: 3,0m >



Welch Allyn

Temperature probe.

Compatible with:

Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.400PT	
Skin TITANIUM® 5,5mm		BIOTEMP.400PTP	
Timpanic		BIOTEMP.400T	
Adult esophageal / rectal		BIOTEMP.400RA	
Pediatric esophageal / rectal		BIOTEMP.400RP	


Welch Allyn

One piece ECG cable (3 or 5 leads)


Compatible with:

Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



< Cable length: 3,5m >




Welch Allyn

ECG trunk cable with leads (3 or 5 leads)


Compatible with:

Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




Welch Allyn

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code
3 leads		BIOECG.BC.3T
5 leads		BIOECG.BC.5T







< Cable length: 3,0m >



Compatible with: Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol.

Welch Allyn


ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code
3 leads snap		BIOECG.RA.S3
3 leads clamp		BIOECG.RA.C3
5 leads snap		BIOECG.RA.S5
5 leads clamp		BIOECG.RA.C5


Compatible with:

Propaq 100/200 Series, Atlas Propaq Encore, Propaq CS, Ultra Smart Protocol.

Cable length: 1,0m >










< Cable length: 3,0m >



YSI 400

Temperature probe.


Options		Code
Skin TITANIUM® 9,5mm		BIOTEMP.400PT
Skin TITANIUM® 5,5mm		BIOTEMP.400PTP
Timpanic		BIOTEMP.400T
Adult esophageal / rectal		BIOTEMP.400RA
Pediatric esophageal / rectal		BIOTEMP.400RP
Air 1		BIOTEMP.400A1
Air 2		BIOTEMP.400A2

Compatible with:

CASMed: 750; Criticare: 506, 507, 507E, 507ELC, 507S, 508, 1100, 2200, 8100 series; Datascope: Passport 2, Spectrum, Spectrum OR, Passport XG, Trio; Datex Ohmeda: S/3, CS/3, S/5, Cardiocap 1, Cardiocap 2, Cardiocap 5, Light Monitor, Satellite Plus; Dräger: Delta, Vista, Gamma, Kappa; Invivo: MDE Prism; Mindray: PM 7000, PM 8000, PM 9000 Express, PM 9000, VS-800; Nihon Kohden: BSM-2300 Life Scope I, BSM-4100 Life Scope P, BSM-5100 Life Scope A, BSM-9510 Life Scope M, BSM-9800 Life Scope S.; Welch Allyn: Atlas.



< Cable length: 3,0m >



YSI 700

Temperature probe.

Options		Code	
Skin TITANIUM® 9,5mm		BIOTEMP.700PT	
Adult esophageal / rectal		BIOTEMP.700RA	
Pediatric esophageal / rectal		BIOTEMP.700RP	


Compatible with: All YSI 700 devices.

Zoll


One piece ECG cable
(3 or 5 leads)

Compatible with:
D900, PD1200, PD1400,
PD1600, PD1700, PD2000.

Options		Code	
3 leads snap		BIOECG.BC.S3F	
3 leads clamp		BIOECG.BC.C3F	
5 leads snap		BIOECG.BC.S5F	
5 leads clamp		BIOECG.BC.C5F	



< Cable length: 3,5m >




Zoll


ECG trunk cable with leads
(3 or 5 leads)

Compatible with:
D900, PD1200, PD1400,
PD1600, PD1700, PD2000.

Options		Code	
3 leads snap		BIOECG.BC.S3D	
3 leads clamp		BIOECG.BC.C3D	
5 leads snap		BIOECG.BC.S5D	
5 leads clamp		BIOECG.BC.C5D	






< Cable length: 4,0m >




Zoll

ECG trunk cable
(3 or 5 leads).
this product requires a lead set

Options		Code	
3 leads		BIOECG.BC.3T	
5 leads		BIOECG.BC.5T	







< Cable length: 3,0m >



Compatible with: D900, PD1200, PD1400, PD1600, PD1700, PD2000.

Zoll

ECG lead set
(3 or 5 leads)
this product requires a trunk cable

Options		Code	
3 leads snap		BIOECG.RA.S3	
3 leads clamp		BIOECG.RA.C3	
5 leads snap		BIOECG.RA.S5	
5 leads clamp		BIOECG.RA.C5	

Compatible with:
D900, PD1200, PD1400, PD1600,
PD1700, PD2000.

Cable length: 1,0m >

