


## AUDIOMETRAI

Eil. Nr.	Aprašymas
1.	 <p><b>SMART 130S</b>  <b>Travel PC Base Tone Screening Audiometer.</b>            Du nepriklausomi kanalai.            Rankinis ir automatinis režimas.            Tone audiometrija –AC-10dB/130dB (PD-81)            SISI TEST 5 ir 1 dB</p>
2.	<p><b>SMART 130P</b>  <b>Travel PC Base Tone Diagnostic Audiometer</b>            Du nepriklausomi kanalai.            Rankinis ir automatinis režimas.            Tone audiometrija –AC-10dB/130dB (PD-81)            BC-10dB/70dB            SISI TEST 5 ir 1 dB</p>



SISI TEST 5 ir 1 dB

5.



**AUDIO 4002**

**PC Base Clinical Audiometer**

Du nepriklausomi kanalai.

Rankinis ir automatinis režimas.

Tone audiometrija –AC-10dB/130dB (PD-81)

BC-10dB/70dB Speech audiometry with world list

SISI TEST 5 ir 1 dB, ABLB, Langebeck