








## EME PHYSIO MAGNETOTERAPIJOS PRIETAISAI

Eil. Nr.	Kodas Modelis	Aprašymas
1.	 <p><b>MA1330 – Magnetomed 2000</b>  <b>MA1330/110 – Magnetomed 2000</b>  <b>110V</b>  <b>(Business line)</b></p>	<p>Equipment for magnetic therapy with 2 independent output channels. Graphic display (240 x 128 pixels) 50 stored protocols 50 new protocols can be stored in the internal memory                      – Frequency 1 - 100 Hz – Power up to 100 Gauss for each channel.</p> <p>Included standard accessories:                      n°1 metal ring to check emission                      n°1 power cable                      n°2 fuses                      n°1 user manual in DVD</p> <p>**For a correct use of the equipment you need to choose one or more accessories from pg. 15</p>
2.	 <p><b>MA1331 – Magnetomed 7200</b>  <b>MA1331/110 – Magnetomed 7200</b>  <b>110V</b>  <b>(Prestige line)</b></p>	<p>Equipment for magnetic therapy with 2 independent output channels. Colour graphic display (320 x 240 pixels) with Touch&amp;Scroll function – 100 stored protocols – internal memory for over 200 new protocols – smart card – frequency 1-100 Hz – Power up to 100 Gauss for each channel.</p> <p>Included standard accessories:                      n°1 metal ring to check emission                      n°1 smart card                      n°1 power cable                      n°2 fuses                      n°1 user manual in DVD</p> <p>**For a correct use of the equipment you need to choose one or more accessories from pg. 15</p>
3.	 <p><b>MA1333 – Magnetomed 8200</b>  <b>MA1333/110 – Magnetomed 8200</b>  <b>110V</b>  <b>(Prestige line)</b></p>	<p>Equipment for magnetic therapy with built-in trolley, 2 independent output channels. Colour graphic display (320 x 240 pixels) with Touch&amp;Scroll function – 100 stored protocols – internal memory for over 200 new protocols – smart card – frequency 1-100 Hz – power up to 100 Gauss for each channel.</p> <p>Included standard accessories:                      n°1 metal ring to check emission                      n°1 smart card                      n°1 power cable                      n°2 fuses                      n°1 user manual in DVD</p> <p>**For a correct use of the equipment you need to choose one or more accessories from pg. 15</p>

4.	 <p><b>MA1334 – Magnetomed 8400</b>  <b>MA1334/110 –Magnetomed 8400</b>  <b>110V</b>  <b>(Prestige line)</b></p>	<p>Equipment for magnetic therapy with built-in trolley, 4 independent output channels.  Colour graphic display (320 x 240 pixels) with Touch&amp;Scroll function – 100 stored protocols – internal memory for over 200 new protocols – smart card – frequency 1-100 Hz – power up to 100 Gauss for each channel.</p> <p>Included standard accessories:  n°1 metal ring to check emission  n°1 smart card  n°1 power cable  n°2 fuses  n°1 user manual in DVD</p> <p><b>**For a correct use of the equipment you need to choose one or more accessories from pg. 15</b></p>
5.	 <p><b>ACC5510 – LC60</b>  <b>Manual sliding</b></p>	<p>High-tech couch in light-weight aluminium with 1 sliding track, 1 cylinder W.I.S (Wood Injection System) <math>\phi</math> 60 cm and 30 cm large.  Manual sliding system</p> <p>Size:  Length 196 cm  Height 70 cm  Width 65 cm</p> <p>Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626).</p>
6.	 <p><b>ACC7150 – LC60</b>  <b>Electric sliding</b></p>	<p>High-tech couch in light-weight aluminium with 1 sliding track, 1 cylinder W.I.S (Wood Injection System) <math>\phi</math> 60 cm and 30 cm large, 1 remote control.  Electric sliding system</p> <p>Size:  Length 196 cm  Height 70 cm  Width 65 cm</p> <p>Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626).</p>
7.	 <p><b>ACC5512 – LC60/2</b></p>	<p>High-tech couch in light-weight aluminium with 2 sliding tracks, 2 cylinders W.I.S (Wood Injection System) <math>\phi</math> 60 cm and 30 cm large. Ideal solution for “total body” treatment.</p> <p>Size:  Length 196 cm  Height 70 cm  Width 65 cm</p> <p>Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626).</p>
8.	ACC5450 – C60	<p>Cylinder W.I.S (Wood Injection System) <math>\phi</math> 60 cm and 30 cm large.  Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626).  Complete with carriage for sliding</p>
9.	ACC5451 – CP60	<p>Portable cylinder W.I.S (Wood Injection System) (with feet) <math>\phi</math> 60 cm.  Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626).</p>
10.	ACC5452 – CP30	<p>Portable cylinder W.I.S (Wood Injection System) (with feet) <math>\phi</math> 30 cm.  Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626).</p>
11.	ACC615 – CP	Pair of applicators for magnetic therapy (16 x 10 x 3,5 cm)
12.	ACC552 – L	High-tech couch in light-weight aluminium

<b>13.</b>	ACC668	Magnetomed cable, linking generator to solenoid
<b>14.</b>	ACC606	Pack of 10 patient smart cards
<b>15.</b>	CONT72	Transport bag in tnt fabric
<b>16.</b>	ACC604	3-shelf trolley (60 x 37 x 86 cm)
		<b>SCIENTIFIC AND MARKETING DOCUMENTATION</b>
<b>17.</b>	ACCM01	“Manuale clinico di terapia fisica” by Dr. Aldo Vasta in Italian (min. 5 copies)
<b>18.</b>	ACCM02	“Clinical Manual Of Physical Therapy” by Dr. Aldo Vasta in English (min. 5 copies)
<b>19.</b>	ACCM04	User manual in paper version
<b>20.</b>	ACCM17	Poster “Magnetomed”