








EME PHYSIO MAGNETOTERAPIJOS PRIETAISAI

| Eil. Nr. | Kodas Modelis | Aprašymas |
|----------|--|---|
| 1. |  <p>MA1330 – Magnetomed 2000 MA1330/110 – Magnetomed 2000 110V (Business line)</p> | <p>Equipment for magnetic therapy with 2 independent output channels. Graphic display (240 x 128 pixels) 50 stored protocols 50 new protocols can be stored in the internal memory – Frequency 1 - 100 Hz – Power up to 100 Gauss for each channel.</p> <p>Included standard accessories: n°1 metal ring to check emission n°1 power cable n°2 fuses n°1 user manual in DVD</p> <p>**For a correct use of the equipment you need to choose one or more accessories from pg. 15</p> |
| 2. |  <p>MA1331 – Magnetomed 7200 MA1331/110 – Magnetomed 7200 110V (Prestige line)</p> | <p>Equipment for magnetic therapy with 2 independent output channels. Colour graphic display (320 x 240 pixels) with Touch&Scroll function – 100 stored protocols – internal memory for over 200 new protocols – smart card – frequency 1-100 Hz – Power up to 100 Gauss for each channel.</p> <p>Included standard accessories: n°1 metal ring to check emission n°1 smart card n°1 power cable n°2 fuses n°1 user manual in DVD</p> <p>**For a correct use of the equipment you need to choose one or more accessories from pg. 15</p> |
| 3. |  <p>MA1333 – Magnetomed 8200 MA1333/110 – Magnetomed 8200 110V (Prestige line)</p> | <p>Equipment for magnetic therapy with built-in trolley, 2 independent output channels. Colour graphic display (320 x 240 pixels) with Touch&Scroll function – 100 stored protocols – internal memory for over 200 new protocols – smart card – frequency 1-100 Hz – power up to 100 Gauss for each channel.</p> <p>Included standard accessories: n°1 metal ring to check emission n°1 smart card n°1 power cable n°2 fuses n°1 user manual in DVD</p> <p>**For a correct use of the equipment you need to choose one or more accessories from pg. 15</p> |

| | | |
|-----------|---|--|
| <p>4.</p> |  <p>MA1334 – Magnetomed 8400 MA1334/110 –Magnetomed 8400 110V (Prestige line)</p> | <p>Equipment for magnetic therapy with built-in trolley, 4 independent output channels. Colour graphic display (320 x 240 pixels) with Touch&Scroll function – 100 stored protocols – internal memory for over 200 new protocols – smart card – frequency 1-100 Hz – power up to 100 Gauss for each channel.</p> <p>Included standard accessories: n°1 metal ring to check emission n°1 smart card n°1 power cable n°2 fuses n°1 user manual in DVD</p> <p>**For a correct use of the equipment you need to choose one or more accessories from pg. 15</p> |
| <p>5.</p> |  <p>ACC5510 – LC60 Manual sliding</p> | <p>High-tech couch in light-weight aluminium with 1 sliding track, 1 cylinder W.I.S (Wood Injection System) ϕ 60 cm and 30 cm large. Manual sliding system</p> <p>Size: Length 196 cm Height 70 cm Width 65 cm</p> <p>Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626).</p> |
| <p>6.</p> |  <p>ACC7150 – LC60 Electric sliding</p> | <p>High-tech couch in light-weight aluminium with 1 sliding track, 1 cylinder W.I.S (Wood Injection System) ϕ 60 cm and 30 cm large, 1 remote control. Electric sliding system</p> <p>Size: Length 196 cm Height 70 cm Width 65 cm</p> <p>Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626).</p> |
| <p>7.</p> |  <p>ACC5512 – LC60/2</p> | <p>High-tech couch in light-weight aluminium with 2 sliding tracks, 2 cylinders W.I.S (Wood Injection System) ϕ 60 cm and 30 cm large. Ideal solution for “total body” treatment.</p> <p>Size: Length 196 cm Height 70 cm Width 65 cm</p> <p>Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626).</p> |

| | | |
|-----|----------------|---|
| 8. | ACC5450 – C60 | Cylinder W.I.S (Wood Injection System) ϕ 60 cm and 30 cm large. Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626). Complete with carriage for sliding |
| 9. | ACC5451 – CP60 | Portable cylinder W.I.S (Wood Injection System) (with feet) ϕ 60 cm. Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626). |
| 10. | ACC5452 – CP30 | Portable cylinder W.I.S (Wood Injection System) (with feet) ϕ 30 cm. Exclusive MFC patented system which concentrates the lines of force inside the solenoid without external dispersion (in compliance with d. lgs. 626). |
| 11. | ACC615 – CP | Pair of applicators for magnetic therapy (16 x 10 x 3,5 cm) |
| 12. | ACC552 – L | High-tech couch in light-weight aluminium |
| 13. | ACC668 | Magnetomed cable, linking generator to solenoid |
| 14. | ACC606 | Pack of 10 patient smart cards |
| 15. | CONT72 | Transport bag in tnt fabric |
| 16. | ACC604 | 3-shelf trolley (60 x 37 x 86 cm) |
| | | SCIENTIFIC AND MARKETING DOCUMENTATION |
| 17. | ACCM01 | “Manuale clinico di terapia fisica” by Dr. Aldo Vasta in Italian (min. 5 copies) |
| 18. | ACCM02 | “Clinical Manual Of Physical Therapy” by Dr. Aldo Vasta in English (min. 5 copies) |
| 19. | ACCM04 | User manual in paper version |
| 20. | ACCM17 | Poster “Magnetomed” |